

## AeroWalk Ultra-Lite Rollator

HCPCS: E0143 + E0156, E0143 + E0156



The AeroWalk series redefines ultra-lightweight rollators, delivering the durability of aluminum with a weight comparable to carbon fiber without the premium price. Designed with a seamless, weld-free frame, it offers a sleek, modern aesthetic while enhancing strength and longevity. Its sturdy side-to-side folding mechanism ensures easy transport, while multiple frame and fork finishes provide a refined, customizable look. Balancing innovation, style, and affordability, the AeroWalk series is the ultimate choice for mobility-focused users seeking premium performance.

- Ultra-Lightweight Design with sleek aluminum frame weighs 14.5 lbs.
- Seamless frame has no visible welds and three sleek color options to choose from.
- Concealed brake cable for improved durability and safety.
- Adjustable handle height 31.5" - 37.3".
- Backrest is adjustable and extends up 37.6" for premium lumbar support.
- Large 10" front wheels provide optimal stability and rolling comfort.
- Durable and comfortable seat for rest on the go.
- Easy side-to-side folding for storage and transport.
- Removable storage bag to easily carry personal items.

**PRODUCTS****Description**

10277BL	Aerowalk Rollator, Blue, 1/ea
10277GY	Aerowalk Rollator, Gray, 1/ea
10277WT	Aerowalk Rollator, White, 1/ea

**SPECIFICATIONS****Description**

Warranty	Limited Lifetime
Product Weight	14.5 lb
Material	Aluminum
Seat Dimensions	17.7" x 10"
Handle Height	31.5" - 37.3"
Weight Capacity	330 lb
Seat To Floor Height	23.1"
Casters	10" Front / 8" Back
Color	Available in: Gray, Blue, White
Dimensions	23.3" x 26.6"
HCPCS	E0143 + E0156
Width Between Arms	19.3"