

Walking Sling

M921503

- Fully adjustable straps
- Simple to apply padded harness
- Additional padding to under arm sections
- Adjustable groin straps with removable comfort covers
- Designed to assist with standing/walking
- Some weight bearing ability recommended
- Wash Temp 95°C/203°F
- Tumble dry low heat (max 50°C)



Product Codes and Specifications:



| 100% Polyester | | A | B | C | D |
|----------------|--------------------|--------------------------|---------------------|--------------|-------------|
| XS | M921503-XS | 30¼ - 34¼" (77 - 87cm) | 8¾" (22cm) or over | 31" (79cm) | 9½" (24cm) |
| S | M921503-S | 34¼ - 42¼" (87 - 107cm) | 10" (25cm) or over | 36¼" (92cm) | 10¾" (27cm) |
| M | M921503-M | 38¼ - 46" (97 - 117cm) | 10¾" (27cm) or over | 40½" (103cm) | 11" (28cm) |
| L | M921503-L | 44 - 52" (112 - 132cm) | 11" (28cm) or over | 45" (114cm) | 11" (28cm) |
| XL | M921503-XL | 49¼ - 57" (125 - 145cm) | 11¾" (30cm) or over | 49½" (126cm) | 11¾" (30cm) |
| XXL | M921503-XXL | 54¼ - 62¼" (138 - 158cm) | 13½" (34cm) or over | 53½" (136cm) | 12½" (32cm) |

