

Company: <input type="text"/>	Ship To: <input type="text"/>
Account Number: <input type="text"/> Date: <input type="text"/>	Address: <input type="text"/>
<input type="checkbox"/> Quote: <input type="checkbox"/> PO Number: <input type="text"/>	City: <input type="text"/> State: <input type="text"/>
ATP: <input type="text"/>	Zip: <input type="text"/> Phone: <input type="text"/>
Marked For: <input type="text"/>	Email/Fax: <input type="text"/>
Associated Quote/Order#: <input type="text"/>	Effective: 1/5/2026

- ☐ Catalyst 5 - Transit Included (K0005) \$2822
- ☐ Catalyst 5 - Transit Not Included (K0005) \$2472
- 300 lb weight capacity

Power Assist

- ☐ Twist Drive System (E0986) \$7950
- Includes Multi-Function Toggle Controller, Removable Li-ion Battery - 36V, 4Ah, Twist Carry Bag, and mounting systems. See page 9 for additional, required selections.*
- 265lb (120 kg) Weight Capacity.

Frame

<h3>Seat Width</h3> <p>21" and 22" not available with Klaxon Twist</p> <table border="0"> <tr><td><input type="checkbox"/> 14"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 15"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 16"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 17"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 18"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 19"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 20" (E2201)</td><td>\$577</td></tr> <tr><td><input type="checkbox"/> 21" (E2201)</td><td>\$577</td></tr> <tr><td><input type="checkbox"/> 22" (E2201)</td><td>\$577</td></tr> </table> <h3>Seat Depth</h3> <table border="0"> <tr><td><input type="checkbox"/> 14"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 15"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 16"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 17"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 18"</td><td>N/C</td></tr> <tr><td><input type="checkbox"/> 19"</td><td>\$81</td></tr> <tr><td><input type="checkbox"/> 20" (E2203)</td><td>\$635</td></tr> </table> <h3>Transit Option</h3> <table border="0"> <tr> <td><input type="checkbox"/> Transit System (Weight Limit 300lbs) (E1022)</td> <td>\$350</td> </tr> <tr> <td colspan="2"><i>Not available with 350 lb Weight Capacity</i></td> </tr> <tr> <td><input type="checkbox"/> Omit</td> <td>N/C</td> </tr> </table>	<input type="checkbox"/> 14"	N/C	<input type="checkbox"/> 15"	N/C	<input type="checkbox"/> 16"	N/C	<input type="checkbox"/> 17"	N/C	<input type="checkbox"/> 18"	N/C	<input type="checkbox"/> 19"	N/C	<input type="checkbox"/> 20" (E2201)	\$577	<input type="checkbox"/> 21" (E2201)	\$577	<input type="checkbox"/> 22" (E2201)	\$577	<input type="checkbox"/> 14"	N/C	<input type="checkbox"/> 15"	N/C	<input type="checkbox"/> 16"	N/C	<input type="checkbox"/> 17"	N/C	<input type="checkbox"/> 18"	N/C	<input type="checkbox"/> 19"	\$81	<input type="checkbox"/> 20" (E2203)	\$635	<input type="checkbox"/> Transit System (Weight Limit 300lbs) (E1022)	\$350	<i>Not available with 350 lb Weight Capacity</i>		<input type="checkbox"/> Omit	N/C	<h3>Weight Capacity</h3> <table border="0"> <tr><td><input type="checkbox"/> 300 lbs</td><td>N/C</td></tr> <tr> <td><input type="checkbox"/> 350 lbs (Heavy Duty)</td> <td>\$588</td> </tr> <tr> <td colspan="2"><i>Requires minimum 5" caster height. Not available with Klaxon Twist.</i></td> </tr> <h3>Cross Brace</h3> <table border="0"> <tr><td><input type="checkbox"/> Standard Cross Braces</td><td>N/C</td></tr> <tr> <td><input type="checkbox"/> Seating System Cross Braces</td> <td>N/C</td> </tr> <tr> <td colspan="2"><i>Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.</i></td> </tr> <tr> <td><input type="checkbox"/> Titanium Cross Braces</td> <td>\$755</td> </tr> <tr> <td colspan="2"><i>Not available with 21" and 22" Seat Widths</i></td> </tr> <h3>Seat Pan</h3> <table border="0"> <tr> <td><input type="checkbox"/> Removable Aluminum Seat Pan E2231</td> <td>\$312</td> </tr> <tr> <td colspan="2"><i>1/2" Shorter than Depth</i></td> </tr> <h3>Rear Frame Type</h3> <p><i>No selection needed with Titanium Frame</i></p> <table border="0"> <tr> <td><input type="checkbox"/> Standard (Curved)</td> <td>N/C</td> </tr> <tr> <td colspan="2"><i>Select from Standard or Angle Adjustable Back Post on chart on the next page</i></td> </tr> <tr> <td><input type="checkbox"/> Depth Adjustable (K0108)</td> <td>\$224</td> </tr> <tr> <td colspan="2"><i>Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.</i></td> </tr> <tr> <td><input type="checkbox"/> Reclining Frame</td> <td>\$1227</td> </tr> <tr> <td colspan="2"><i>Select from Reclining Back Post in chart on the next page. Not available with Dynamic Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.</i></td> </tr> </table> </table></table></table>	<input type="checkbox"/> 300 lbs	N/C	<input type="checkbox"/> 350 lbs (Heavy Duty)	\$588	<i>Requires minimum 5" caster height. Not available with Klaxon Twist.</i>		<input type="checkbox"/> Standard Cross Braces	N/C	<input type="checkbox"/> Seating System Cross Braces	N/C	<i>Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.</i>		<input type="checkbox"/> Titanium Cross Braces	\$755	<i>Not available with 21" and 22" Seat Widths</i>		<input type="checkbox"/> Removable Aluminum Seat Pan E2231	\$312	<i>1/2" Shorter than Depth</i>		<input type="checkbox"/> Standard (Curved)	N/C	<i>Select from Standard or Angle Adjustable Back Post on chart on the next page</i>		<input type="checkbox"/> Depth Adjustable (K0108)	\$224	<i>Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.</i>		<input type="checkbox"/> Reclining Frame	\$1227	<i>Select from Reclining Back Post in chart on the next page. Not available with Dynamic Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.</i>	
<input type="checkbox"/> 14"	N/C																																																																						
<input type="checkbox"/> 15"	N/C																																																																						
<input type="checkbox"/> 16"	N/C																																																																						
<input type="checkbox"/> 17"	N/C																																																																						
<input type="checkbox"/> 18"	N/C																																																																						
<input type="checkbox"/> 19"	N/C																																																																						
<input type="checkbox"/> 20" (E2201)	\$577																																																																						
<input type="checkbox"/> 21" (E2201)	\$577																																																																						
<input type="checkbox"/> 22" (E2201)	\$577																																																																						
<input type="checkbox"/> 14"	N/C																																																																						
<input type="checkbox"/> 15"	N/C																																																																						
<input type="checkbox"/> 16"	N/C																																																																						
<input type="checkbox"/> 17"	N/C																																																																						
<input type="checkbox"/> 18"	N/C																																																																						
<input type="checkbox"/> 19"	\$81																																																																						
<input type="checkbox"/> 20" (E2203)	\$635																																																																						
<input type="checkbox"/> Transit System (Weight Limit 300lbs) (E1022)	\$350																																																																						
<i>Not available with 350 lb Weight Capacity</i>																																																																							
<input type="checkbox"/> Omit	N/C																																																																						
<input type="checkbox"/> 300 lbs	N/C																																																																						
<input type="checkbox"/> 350 lbs (Heavy Duty)	\$588																																																																						
<i>Requires minimum 5" caster height. Not available with Klaxon Twist.</i>																																																																							
<input type="checkbox"/> Standard Cross Braces	N/C																																																																						
<input type="checkbox"/> Seating System Cross Braces	N/C																																																																						
<i>Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.</i>																																																																							
<input type="checkbox"/> Titanium Cross Braces	\$755																																																																						
<i>Not available with 21" and 22" Seat Widths</i>																																																																							
<input type="checkbox"/> Removable Aluminum Seat Pan E2231	\$312																																																																						
<i>1/2" Shorter than Depth</i>																																																																							
<input type="checkbox"/> Standard (Curved)	N/C																																																																						
<i>Select from Standard or Angle Adjustable Back Post on chart on the next page</i>																																																																							
<input type="checkbox"/> Depth Adjustable (K0108)	\$224																																																																						
<i>Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.</i>																																																																							
<input type="checkbox"/> Reclining Frame	\$1227																																																																						
<i>Select from Reclining Back Post in chart on the next page. Not available with Dynamic Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.</i>																																																																							

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0011.47 Pricing and specifications subject to change without notice.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 6). Requires Rear Anti-Tips. Not available with Spinergy wheels or Height Adjustable Flip Back T-Arm.

☐ e-motion Ready N/C

Titanium Option

☐ Titanium Frame \$2073
Includes hangers and base frame. Does not include Cross Braces. Not available with Frame Color.

Backrest

Depth Adjustable not available with Titanium Frame

Back Post		Price	Height Range	Back Height				
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> Short ^{3,4}	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10.5"	<input type="checkbox"/> 11.5"		
		<input type="checkbox"/> Medium ⁴	<input type="checkbox"/> 12.5" ⁶	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15.5"		
		<input type="checkbox"/> Tall	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19.5"		
<input type="checkbox"/> Angle Adjustable w/ Push Handle	\$208	<input type="checkbox"/> Short ⁴	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"			
		<input type="checkbox"/> Medium ⁴	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"			
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"			
<input type="checkbox"/> Stroller Handle Back Post (Depth Adjustable)	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 20"					
		<input type="checkbox"/> Tall	<input type="checkbox"/> 24"					
<input type="checkbox"/> Height Adjustable Straight with Push Handle (Depth Adjustable) ^{2,9}	\$101	<input type="checkbox"/> Short ¹	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	
		<input type="checkbox"/> Tall	<input type="checkbox"/> 20" ⁸	<input type="checkbox"/> 21" ⁵	<input type="checkbox"/> 22" ⁵	<input type="checkbox"/> 23" ⁵	<input type="checkbox"/> 24" ⁵	
<input type="checkbox"/> Reclining w/ Adjustable Handle (E1225) ^{1,7}	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 22"					
		<input type="checkbox"/> Tall	<input type="checkbox"/> 26"					

¹ Not available with Dynamic Back

² Not available with Standard Back Upholstery

³ Not available with Tension Adjustable Back Upholstery

⁴ Not available with Standard Back Upholstery in combination with 21" and 22" Seat Widths

⁵ Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

⁶ Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

⁷ Not available with O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

⁸ Not available with Angle Adjustable Locking Extendable Flip Up Armrest

⁹ Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrests

Reclining Back Options

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

☐ Standard Spring N/C
☐ Heavy Duty Spring N/C

Standard Back Post

Select option below only if a Standard Back Post option is selected

☐ Straight with Push Handle N/C
☐ Straight without Push Handle N/C
☐ 8 Degree Bend with Push Handle \$37
Not available with Short Back Height Range
☐ Straight with Fold Down Push Handle \$134
Not available with 350 lb Weight Capacity
☐ Caregiver Adjustable Height \$167
Not available with Swing Away Armrests. Available only with Straight without Push Handle Backrest.

Push Handle Options

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 lb Weight Capacity

☐ Low Profile Push Handle \$96
2" shorter than standard push handles

Depth Adjustable Back Post Options

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

☐ Stroller Handle Extensions \$101
Not available with Height Adjustable Straight with Push Handle Back Post
☐ Dynamic Back (E2398) \$1451
Available only with T-Arms

Back Upholstery Options

Not available with Reclining or Stroller Handle Back Post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

☐ Standard N/C
☐ Tension Adjustable with Standard Cover (E2611) \$359
☐ Tension Adjustable with Air Mesh Breathable Cover (E2611) \$422
☐ Omit (\$15)

U.S. Military Service Patches

Attached only if ordered with Tension Adjustable Back Upholstery

☐ Army ☐ Navy ☐ Marines ☐ Air Force
☐ Coast Guard ☐ US Purple Heart \$96
☐ Combo Army/Purple Heart ☐ Combo Navy/Purple Heart
☐ Combo Marines/Purple Heart ☐ Combo Air Force/Purple Heart
☐ Combo Coast Guard/Purple Heart \$155

Frame Type / Footrest

Frame Type	
<input type="checkbox"/> Swing Away	N/C
<input type="checkbox"/> Super Low - Swing Away	N/C
<input type="checkbox"/> 70° Fixed Front	N/C
<input type="checkbox"/> 80° Fixed Front	N/C

ELR Pad and Pro ELR Cover	
<input type="checkbox"/> Standard ELR Foam Knee Pad (Pair) K0108	\$386
<i>Only available with ELR</i>	
<input type="checkbox"/> Pro ELR Gel Pad Cover (Pair) K0108	\$468
<i>Only available with Pro ELR</i>	

Footrest Hangers

Frame	Hanger	60°	70°	80°	90°		Retail
Swing Away	Extension Tube Mount	<input type="checkbox"/> ²	<input type="checkbox"/>	<input type="checkbox"/> ¹			N/C
	Front Mount ²		<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/> ¹		N/C
	Elevating Leg Rests ^{3,4}					<input type="checkbox"/>	\$307
	Pro Elevating Leg Rests E0990 ^{3,4}					<input type="checkbox"/>	\$354
	Omit Hangers					<input type="checkbox"/>	(\$20)
Fixed	Extension Tube Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹			N/C

¹ Not available with 7" or 8" casters² Not available with Titanium Frame or Klaxon Twist³ Available only with Classic Hanger Release⁴ Not available in combination with 24" wheels and Push to Lock Flush Mount Wheel Locks

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°,70°,80°)	<input type="checkbox"/> Short ³	<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"			
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"			
Front Mount ^{2,3} Hangers		Standard Selections									Requires 3" Extension		
	<input type="checkbox"/> Short Front Mount Hanger (70°, 80°, 90°)	<input type="checkbox"/> 4" ¹	<input type="checkbox"/> 4.5" ¹	<input type="checkbox"/> 5"	<input type="checkbox"/> 5.5"	<input type="checkbox"/> 6"	<input type="checkbox"/> 6.5"				<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"	<input type="checkbox"/> 11"
		<input type="checkbox"/> 7"	<input type="checkbox"/> 7.5"	<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 9.5"				<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"
	<input type="checkbox"/> Medium Front Mount Hanger (70°, 80°, 90°)	<input type="checkbox"/> 4" ¹	<input type="checkbox"/> 4.5" ¹	<input type="checkbox"/> 5"	<input type="checkbox"/> 5.5"	<input type="checkbox"/> 6"	<input type="checkbox"/> 6.5"	<input type="checkbox"/> 7"	<input type="checkbox"/> 7.5"		<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"
		<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"	<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"		<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"
	<input type="checkbox"/> Long Front Mount Hanger (70°, 80°)	<input type="checkbox"/> 5"	<input type="checkbox"/> 5.5"	<input type="checkbox"/> 6"	<input type="checkbox"/> 6.5"	<input type="checkbox"/> 7"	<input type="checkbox"/> 7.5"	<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"
		<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"	<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"
ELR ³	<input type="checkbox"/> Standard ELR Extension Tube Short (11.75" - 16.75")												
	<input type="checkbox"/> Standard ELR Extension Tube Medium (13.75" - 18.75")												
	<input type="checkbox"/> Standard ELR Extension Tube Long (16.25" - 21.25")												
Pro ELR Low ³ Pivot Setting	<input type="checkbox"/> Pro ELR Low Pivot Setting Short (12" - 18")												
	<input type="checkbox"/> Pro ELR Low Pivot Setting Long (14.5" - 21.5")												
Pro ELR High ³ Pivot Setting	<input type="checkbox"/> Pro ELR High Pivot Setting Short (10.5" - 16")												
	<input type="checkbox"/> Pro ELR High Pivot Setting Long (13" - 19.5")												
Fixed Front Frame (70°,80°)	<input type="checkbox"/> Short ³	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"				
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"			

¹ Available with 90° only² See Minimum Front Seat Height chart below³ Not available with Klaxon Twist

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

3" Extension

Available only with Front Mount Hangers

<input type="checkbox"/> 3" Extension	\$38
<input type="checkbox"/> 3" Extension (Provided as Parts)	\$38

Footplates

<input type="checkbox"/> Composite	N/C
<i>Not available with Front Mount Hanger</i>	
<input type="checkbox"/> Composite Angle Adjustable (K0040)	\$208
<i>Not available with Front Mount Hanger or Seat Widths less than 16"</i>	
<input type="checkbox"/> Aluminum Angle Adjustable (K0040)	\$214
<input type="checkbox"/> Aluminum Locking Multi-Angle Adjustable Footplate (K0108)	\$337
<i>Available only with Front Mount Hanger</i>	
<input type="checkbox"/> One Piece Flip Up Angle Adjustable (K0040)	\$208
<i>Not available with Front Mount Hanger, ELR, Pro ELR, or Klaxon Twist</i>	

Hanger Release

If no selection is made Classic Hanger Release will be used

<input type="checkbox"/> Classic	N/C
<input type="checkbox"/> 4-Way - In-Line	N/C
<input type="checkbox"/> 4-Way - Out	N/C

Calf Straps

<input type="checkbox"/> Velcro Adjustable (K0038)	\$43
<input type="checkbox"/> Padded Velcro Adjustable (K0038)	\$75
<input type="checkbox"/> Padded Airmesh Calf Strap (K0038)	\$139

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

<input type="checkbox"/> Heel Loops (E0951)	\$44
<input type="checkbox"/> Adjustable Heel Loops (E0951)	\$54
<input type="checkbox"/> Omit Heel Loops	N/C

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size				
Left	<input type="checkbox"/> Left Limb Support (E1020)	\$449	<input type="checkbox"/> Omit Left Hanger	<input type="checkbox"/> Smooth Cover Left	<input type="checkbox"/> Air Mesh Cover Left	<input type="checkbox"/> 5x8 Pad Left Side	<input type="checkbox"/> 7x10 Pad Left Side	<input type="checkbox"/> 9x10 Pad Left Side	<input type="checkbox"/> 9x14 Pad Left Side	<input type="checkbox"/> 11x10 Pad Left Side
Right	<input type="checkbox"/> Right Limb Support (E1020)	\$449	<input type="checkbox"/> Omit Right Hanger	<input type="checkbox"/> Smooth Cover Right	<input type="checkbox"/> Air Mesh Cover Right	<input type="checkbox"/> 5x8 Pad Right Side	<input type="checkbox"/> 7x10 Pad Right Side	<input type="checkbox"/> 9x10 Pad Right Side	<input type="checkbox"/> 9x14 Pad Right Side	<input type="checkbox"/> 11x10 Pad Right Side

Front Seat Height and Casters

Front Seat Height

Select Front Seat Height in row corresponding to desired caster size

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"
5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	
6"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	<input type="checkbox"/> 21"		
7"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"		<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	<input type="checkbox"/> 21"				
8"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	<input type="checkbox"/> 21"						

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Casters

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C	<input type="checkbox"/>				
	.75" Lighted Rollerblade	\$81	<input type="checkbox"/>				
	1.25" Soft Roll	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1" Poly Aluminum	\$123	<input type="checkbox"/>	<input type="checkbox"/>			
	1" Pneumatic	\$96			<input type="checkbox"/>		<input type="checkbox"/>
	1.5" Soft Roll	\$63	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	1.5" Soft Roll Aluminum	\$221	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	2" Poly	\$101			<input type="checkbox"/>		
	2" Pneumatic w/Foam Insert	\$256					<input type="checkbox"/>

Caster Housing

- ☐ Standard Housing N/C
- ☐ Performance Housing \$54
Not available with 350 lb Weight Capacity

Caster Pin Locks

- ☐ Caster Pin Locks (K0073) \$219
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock

- ☐ Push to Lock N/C
- ☐ Push to Lock Flush Mount N/C
- ☐ Pull to Lock N/C
- ☐ Low Profile Scissor Lock \$75
Not available with low poly or full poly tires
- ☐ Drum Brake (E2228) \$619
Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.
- ☐ Dual - Drum Brake with Push to Lock (E2228) \$731
Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.
- ☐ Dual - Drum Brake with Pull to Lock (E2228) \$731
Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.

Fork Options

- ☐ Best Fit N/C
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
- ☐ Standard Fork Stem N/C
- ☐ 3/4" Longer Fork Stem N/C
- ☐ 1 1/2" Longer Fork Stem N/C
- ☐ Frog Legs Suspension Fork \$508
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic with Foam Insert caster.
- ☐ Frog Legs II Suspension System \$508
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic with Foam Insert caster.

Frog Legs Options

- ☐ User Weight Under 80 lbs N/C
- ☐ User Weight 80-120 lbs N/C
- ☐ User Weight 120-160 lbs N/C
- ☐ User Weight 160-220 lbs N/C
- ☐ User Weight 220-260 lbs N/C
- ☐ User Weight 260-300 lbs N/C

Wheel Lock Options

- ☐ Hemi Wheel Lock - Left (K0108) \$374
Extension Handle included
- ☐ Hemi Wheel Lock - Right (K0108) \$374
Extension Handle included
- ☐ Omit N/C
- ☐ Extension Handles (E0961) \$85
Available only with Push to Lock or Pull to Lock Wheel Locks
- ☐ Grade Aids (E0974) \$166
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires
- ☐ Automatic Braking System (E2228) \$2060
PRATIKO Automatic Braking System. Installed by Ki. Wheel Lock selection required - If no selection, Push to Lock will be provided. Weight capacity of 250 lbs (113 kg). Not available with transit, drum brake, seat pan, OAD, Reclining Back, Heavy Duty, Klaxon Twist, or seat widths wider than 21". Available with the following rear wheels: 24" Spoke, Maxx Spoke and Mag. 22" is Mag only.

Rear Seat Height and Wheels

Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height									
20" ⁴	Low	<input type="checkbox"/> 11.5" ³	<input type="checkbox"/> 12" ¹	<input type="checkbox"/> 12.5" ³	<input type="checkbox"/> 13" ¹	<input type="checkbox"/> 13.5" ³	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5" ³	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5" ³	<input type="checkbox"/> 16" ³
	Full	<input type="checkbox"/> 12" ³	<input type="checkbox"/> 12.5" ¹	<input type="checkbox"/> 13" ³	<input type="checkbox"/> 13.5" ¹	<input type="checkbox"/> 14" ³	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15" ³	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16" ³	<input type="checkbox"/> 16.5" ³
22" ⁵	Low	<input type="checkbox"/> 12.5" ³	<input type="checkbox"/> 13" ¹	<input type="checkbox"/> 13.5" ³	<input type="checkbox"/> 14" ¹	<input type="checkbox"/> 14.5" ³	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5" ³	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5" ³	<input type="checkbox"/> 17" ³
	Full	<input type="checkbox"/> 13" ³	<input type="checkbox"/> 13.5" ¹	<input type="checkbox"/> 14" ³	<input type="checkbox"/> 14.5" ¹	<input type="checkbox"/> 15" ³	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16" ³	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17" ³	<input type="checkbox"/> 17.5" ³
24" ⁶	Low	<input type="checkbox"/> 13.5" ³	<input type="checkbox"/> 14" ¹	<input type="checkbox"/> 14.5" ³	<input type="checkbox"/> 15" ¹	<input type="checkbox"/> 15.5" ³	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5" ³	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5" ³	<input type="checkbox"/> 18" ³
	Full	<input type="checkbox"/> 14" ³	<input type="checkbox"/> 14.5" ¹	<input type="checkbox"/> 15" ³	<input type="checkbox"/> 15.5" ¹	<input type="checkbox"/> 16" ³	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17" ³	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18" ³	<input type="checkbox"/> 18.5" ³
25" ⁷	Low	<input type="checkbox"/> 14" ³	<input type="checkbox"/> 14.5" ¹	<input type="checkbox"/> 15" ³	<input type="checkbox"/> 15.5" ¹	<input type="checkbox"/> 16" ³	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17" ³	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18" ³	<input type="checkbox"/> 18.5" ³
26" ⁸	Low	<input type="checkbox"/>	<input type="checkbox"/> 15" ^{1,2}	<input type="checkbox"/> 15.5" ³	<input type="checkbox"/> 16" ¹	<input type="checkbox"/> 16.5" ³	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5" ³	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5" ³	<input type="checkbox"/> 19" ³
	Full	<input type="checkbox"/> 15" ^{2,3}	<input type="checkbox"/> 15.5" ^{1,2}	<input type="checkbox"/> 16" ³	<input type="checkbox"/> 16.5" ¹	<input type="checkbox"/> 17" ³	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18" ³	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19" ³	<input type="checkbox"/> 19.5" ³

1 Not available with Amputee Axle Plate

2 Not available with Height Adjustable Flip Back T-Arm

3 Requires Universal Axle Plate

4 Not available with Klaxon Twist

5 Lowest Rear Seat Height with Klaxon Twist and 22" wheels is 15" with Low Profile tires and 15.5" with Full Profile tires

6 Lowest Rear Seat Height with Klaxon Twist and 24" wheels is 16" and 16.5" with Mountain tire

7 Lowest Rear Seat Height with Klaxon Twist and 25" wheels is 16" and 17" with Mountain tire

8 Lowest Rear Seat Height with Klaxon Twist and 26" wheels is 17"

One Arm Drive

- ☐ One Arm Drive - Quick Release (E0958) \$1168
Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.
Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber or with Klaxon Twist.

Drive Side

- ☐ Left Side Drive N/C
☐ Right Side Drive N/C

Drive Handrim

- ☐ Aluminum Drive Handrim N/C
☐ Poly Coated Drive Handrim \$219

Axle Plates

No selection needed with One Arm Drive

- ☐ Universal N/C
☐ Amputee (E0959) \$195
Not available with Swing Away Armrests or Wheelchair Drive Ready Option.
Required with Reclining Backrest.

Handrims

- ☐ Aluminum Anodized N/C
☐ Superlight Aluminum Anodized \$59
Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels. Not available with One Arm Drive.
☐ Plastic Coated \$115
☐ Projection (8 Vertical Projections) \$374
Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used

- ☐ Close Mount \$33
Reduces overall width by .75"

Rear Wheel Type

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Maxx Spoke	\$224	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Maxx Spoke	N/C			<input type="checkbox"/>		
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Colored Mag ³	\$125	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Superlight Spoke ¹	\$294			<input type="checkbox"/>	<input type="checkbox"/> ²	
	Spinergy Spox	\$1014		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Spinergy LX	\$1273		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

1 Not available with 350 lb Weight Capacity

2 Available only with Superlight and Plastic Coated Handrims

3 Color selection on page 8

Spinergy Options

Select spoke color for Spinergy Wheels

- ☐ Black Spokes ☐ Red Spokes ☐ Blue Spokes ☐ White Spokes N/C

Ergonomic Handrims

Not available with 20" or 26" wheels

- ☐ Lite Regular Grip \$497
☐ Standard Regular Grip \$475
☐ Lite Extra Grip \$550
☐ Standard Extra Grip \$528
☐ Lite Omit Grip \$422
☐ Standard Omit Grip \$369

Rear Tire Type

Full Profile

- ☐ Street Pneumatic 1 3/8" (E2211, E2212) N/C
Not available on 25" wheel size
- ☐ Full Poly 1 3/8" (E2218) N/C
Not available on Superlight wheel or 25" or 26" wheel sizes
- ☐ Pneumatic with Airless Insert 1 3/8" (E2211, E2213) \$97
Not available on Superlight wheel or 25" wheel size
- ☐ Schwalbe HS-127 1 3/8" (E2211, E2212) \$139
Available only on 24" wheel size
- ☐ Schwalbe Downtown 1 3/8" (E2211, E2212) \$108

Low Profile

- ☐ Low Poly 1" (E2218) N/C
Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Capacity
- ☐ Sentinel Solid Tire 1" (E2218) \$108
Not available on Superlight Wheels
- ☐ High Pressure 1" (E2211, E2212) \$108
Not available on Mag wheel
- ☐ Schwalbe One 1" (E2211, E2212) \$150
Not available on Mag wheel, 20", 22" or 26" wheel sizes
- ☐ Schwalbe Marathon Plus Evolution 1" (E2211, E2212) \$161
Not available on Mag wheel, 20" or 22" wheel sizes

Axles

No selection needed with One Arm Drive or Drum Brake

- ☐ Quick Release N/C
- ☐ Quad Quick Release \$193
Not available with Mag or Spinergy wheels

Axle Position

Rear Wheel Camber

2° and 4° Not available with Heavy Duty option. 4° Not available with Klaxon Twist.

- ☐ 0° ☐ 2° ☐ 4° N/C

Center of Gravity

No selection needed with Amputee Axle Plate

- ☐ Rear ☐ Center ☐ Forward N/C

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk <input type="checkbox"/>	\$337
	Full <input type="checkbox"/>	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk <input type="checkbox"/>	\$337
	Full <input type="checkbox"/>	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk <input type="checkbox"/>	\$337
	Full <input type="checkbox"/>	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973 ²	Desk <input type="checkbox"/>	\$352
	Full <input type="checkbox"/>	
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973 ²	Desk <input type="checkbox"/>	\$352
	Full <input type="checkbox"/>	

¹ If requested configuration is unavailable, Standard Height T-Arm will be used

² Available only with Curved Rear Frame

³ Not available with Seat Depths less than 16"

⁴ Available only with Standard Back Posts

⁵ Available only with Depth Adjustable and Reclining Back Posts

⁶ Not available with Arm Pads

⁷ Not Available with 15" Rear Seat Height and 24" wheels

Armrest Type	Length	Price
Flip Back (10.5") ^{3,4,7}	Desk <input type="checkbox"/>	N/C
	Full <input type="checkbox"/>	
Angle Adjustable Locking Extendable Flip Up K0020 ⁵	Desk <input type="checkbox"/>	\$304
	Full <input type="checkbox"/>	
Swing Away (10.5"-12.5") E0973 ^{4,6}	<input type="checkbox"/>	\$251
Swing Away with Side Guard (10.5"-12.5") E0973 ^{4,6}	<input type="checkbox"/>	\$374
Omit Armrests	<input type="checkbox"/>	N/C

Arm Pads

- ☐ Standard Arm Pad N/C
- ☐ Waterfall Arm Pad \$28
Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low
- ☐ Foam Grip N/C
Available only with Angle Adjustable Locking Extendable Flip Up Armrest

Accessories

Removable Side Guard

Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests

- ☐ Adult Composite Side Guards (K0108) \$161
- ☐ Pediatric Composite Side Guards (K0108) \$161
- ☐ Aluminum Straight - Black Matte (K0108) \$219
- ☐ Aluminum Straight - Match Frame Color (K0108) \$272
Not available with Hydrographics or Titanium Frame
- ☐ Aluminum Adjustable Fender w/Foam Pad - Black Matte (K0108) \$304
- ☐ Aluminum Adjustable Fender w/Foam Pad - Match Frame Color (K0108) \$357
Not available with Hydrographics or Titanium Frame
- ☐ Carbon Fiber Adjustable Fender (K0108) \$534

Anti-tips

- ☐ Rear Anti-Tips (E0971) \$172
- ☐ Swing-In Anti-Tip Single - Left (E0971) \$171
Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.
- ☐ Swing-In Anti-Tip Single - Right (E0971) \$171
Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.
- ☐ Swing-In Anti-Tip Pair (E0971) \$342
Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with O2 Holder.
- ☐ Omit N/C

Pelvic Positioning Belts

<input type="checkbox"/>	1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$91
<input type="checkbox"/>	1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$183
<input type="checkbox"/>	2" Airline Buckle Pelvic Positioning Belt (E0978)	\$120
<input type="checkbox"/>	1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$183
<input type="checkbox"/>	Small Flex-Mount Positioning Belt (16-24") (E0978)	\$208
<input type="checkbox"/>	Medium Flex-Mount Positioning Belt (19-29") (E0978)	\$208
<input type="checkbox"/>	Large Flex-Mount Positioning Belt (25-37") (E0978)	\$208
<input type="checkbox"/>	Belt Mounting Band Clamp	\$37

Clamp includes hardware to mount Belt

O2 Holder

Not Available with Depth Adjustable Rear Frame

<input type="checkbox"/>	O2 Holder - Left Side (E2208)	\$325
<input type="checkbox"/>	O2 Holder - Right Side (E2208)	\$325

Additional Accessories

<input type="checkbox"/>	Tool Kit	\$63
<input type="checkbox"/>	Removable Underseat Pouch	\$64
<input type="checkbox"/>	Pediatric Ki Mobility Backpack	\$74
<input type="checkbox"/>	Ki Mobility Backpack	\$83
<input type="checkbox"/>	Luggage Carrier	\$90

Not available with ELR, Pro ELR or Residual Limb Support

<input type="checkbox"/>	Black Neoprene Impact Guards	\$99
<input type="checkbox"/>	Clear Spoke Guards (K0065)	\$125
<input type="checkbox"/>	Black Spoke Guards (K0065)	\$128
<input type="checkbox"/>	Cane and Crutch Holder	\$214
<input type="checkbox"/>	IV Pole	\$303

Available only with Standard Straight Back Posts

Upper Extremity Support Tray

3/8" Polycarbonate with ABS ribs. Hardware ordered separately. BC = Body Cut Out. Armrest assembly is required for tray mounting.

Tray	Fits Chair	Retail
<input type="checkbox"/> Tray, Small 22" x 20" (BC 13" x 8") (E0950)	10" - 14"	\$235
<input type="checkbox"/> Tray, Medium 24" x 22" (BC 15" x 9") (E0950)	12" - 16"	\$235
<input type="checkbox"/> Tray, Large 26" x 24" (BC 17" x 10") (E0950)	14" - 18"	\$235
<input type="checkbox"/> Tray, X-Large 28" x 26" (BC 18" x 10") (E0950)	16" - 20"	\$255
<input type="checkbox"/> Tray Hardware Assy, Knob Release - Pair		\$222

Color

Touch-up paint is available for purchase through the parts manual

Options	Price	Colors																			
		Gloss														Matte					
		Black	Ruby Red	Radiant Orange	Yellow	Electric Green	Rich Teal	Glacier Blue	Blue Steel	Deep Sapphire Blue	Purple Plum	Opalescent Violet	Dark Pink	Light Pink	Champagne	Shadow Black	Iridescent Navy	Bright Red	Tactical Green	Dry Earth Tan	Platinum
Frame	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Mag Color		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Color Accent Package ¹	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Color Anodizing Package ²	Black	Blue	BrightGreen	Orange	Pewter	Pink	Purple	Red	Gold												
	N/C	\$255	\$255	\$255	\$255	\$255	\$255	\$255	\$255												
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												

¹ Includes axle plates and forks. When One Arm Drive or Amputee Axle Plate is selected, only forks will be painted. When Frog Legs Suspension Fork or Frog Legs II Suspension System is selected, only axle plate will be painted.

² Select desired color. If no selection is made, black will be used. Includes standard caster housing, caster fork, and axle plates. Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System. Not available with amputee axle plate. Extended lead time may apply.

Hydrographics Options

Not available with Frame Color. Hydrographics Base Color required with Hydrographics. Zebra not available with black.

Options	Price	Hydrographics						Hydrographics Base Color					
		AOR1 (Desert Digital Camo)	NEXT Camo	Pink Camo	Carbon Fiber	Muddy Girl	Zebra	Black	Blue	Granny Smith	Pink	Purple	Red
Frame	\$859	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
Base Color	N/C							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Power Assist

The following options are specific to the Twist Power Assist Add-On option

Controller Mounting Position

If no selection is made Right Side Position will be used

- ☐ Left Side N/C
- ☐ Right Side N/C

Folding Hardware Option

Required selection for folding chair installation

- ☐ Folding Hardware Option \$295

Additional Accessories

- ☐ Attendant Controller (K0108) \$315
- ☐ Handlebar Controller (K0108) \$495
- ☐ Twist Carry Bag \$195
- ☐ Handlebar and Front Mount Connector Carry Bag \$195
- ☐ Removable Li-ion Battery - 36V, 4Ah \$595

Installation Options

- ☐ Install Both Front and Rear Hardware N/C
- ☐ Install Front Hardware Only N/C
- ☐ Install Rear Hardware Only N/C
- ☐ Not Installed N/C

Setup Options

For Custom Settings complete the chart(s) below. Selection made indicates maximum performance level. Selection should be made according to client capabilities, use environment, etc. Setting can be changed in Klaxon app. If no selection is made Twist will be set to defaults presented below.

- ☐ Default Settings N/C
- ☐ Custom Settings N/C

Acceleration Rate										
Defines how quickly the Twist accelerates or decelerates in Variable Speed Mode. Lower values are ideal for indoor or tight spaces. Higher values offer faster acceleration, better suited for outdoor use.										
Applicable Modes: Front and Rear										
Default: 2										
<input type="checkbox"/> 0	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
Maximum Speed										
Sets the top speed limit for the Twist in Variable Speed Mode. If front-mounted, a safe speed limit is set at 4 mph. Rear-mounted increases it to 6 mph.										
Applicable Modes: Front and Rear										
Default: 6 mph										
	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6				
Manual Rim Braking - Force Required										
Adjusts the force required to manually brake using the wheelchair's rims. Lower values mean less force is needed to initiate braking, but may result in false detections over bumps or gaps. Higher values require more force to trigger braking.										
Applicable Modes: Front and Rear										
Default: 5										
<input type="checkbox"/> 0	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
Turning Speed Reduction										
Adjusts how much the Twist reduces speed when steering.										
Applicable Modes: Front										
Default: 0										
<input type="checkbox"/> 0	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
Adaptive Drive Control										
Adjusts how the Twist dynamically responds to user weight, input and terrain. Lower values improve control; higher provide more power.										
Applicable Modes: Front and Rear										
Default: 1										
	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10

Additional custom settings are available in KiO and on the Twist Order form, located at www.kimobility.com

Catalyst Front Seat Height Matrix																					
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster					
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork		
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"		
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4			
15.5 ¹	14 ¹	17.5		16	1 1/2	18.5		17	1 1/2	18.5		17.5	20.5	n/a		1 1/2					
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	3/4		19	17.5	Std	7"		
16.5 ¹	15 ¹	1 1/2		16.5	15			17.5	16			Std	19.5			18	19.5			18	
17	15.5	3/4		16.5	15.5			18	16.5			20	18.5			20	18.5				
15.5	14	Std		17	16	1 1/2		17.5	16.5	3/4		6"	18.5	17		Std	7"			20.5	19
16	14.5			17.5	16			18	17				19	17.5						19.5	18.5
16	15	3/4		18	16.5	18.5		17.5	19.5				18	20						19	
16.5	15.5	1 1/2	16.5	15	Std	6"		18.5	17	1 1/2			7"	20	18.5	3/4		20.5		n/a	3/4
17	15.5		17	15.5				19	17.5					19	18			21		n/a	
17.5	16	17.5	16	19.5	18			19.5	18.5	20.5				19	1 1/2						
16.5	15	Std	17.5	16	3/4		6"	18	16.5	Std	7"			20	19	3/4		21		n/a	1 1/2
17	15.5		18	16.5				18.5	17					20.5	n/a						
17	16	3/4	18.5	17	1 1/2			19	17.5	1 1/2				20	18.5						
17.5	16.5	18	16.5	19.5				18	20.5			19									
18	16.5	18.5	17	18.5				17.5	21			n/a									
18.5	17	1 1/2	19	17.5	3/4			7"	19	18		3/4		7"							
18	16.5		18	16.5		19.5			18.5												
18.5	17	18.5	17	20	19	1 1/2															
18.5	17.5	19	17.5	19.5	18																
19	18	18.5	17.5	20	18.5																
19.5	18.5	1 1/2	19	18	3/4	20.5	19		1 1/2												
20	19		19.5	18.5		21	n/a														
				19.5	18	1 1/2	6" x 2" caster not available with 5" fork														
				20	18.5																
				20.5	19																

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.