

# JEFFERSON SPECIFICATIONS

As always, we accept returns and  
charge NO restocking fees!

See DEMO VIDEOS:  
[www.youtube.com/user/MerrittManufacturing](https://www.youtube.com/user/MerrittManufacturing)

Minimum weight: 7.5 lbs.  
Maximum weight: 40 lbs.

Each Jefferson car seat comes with two  
support cushions, one for each buckle po-  
sition, to optimally position smaller children.

Two buckle positions ensure buckle/  
harness does not come into contact with  
omphalocele.

Soft foam buckle pad to protect child from  
the buckle.

Durable medical-grade naugahyde  
cleanable cover.



Designed and manufactured in the United States by  
Merritt Manufacturing, Inc. The Jefferson™ conforms to all  
applicable Federal motor vehicle safety standards.

©2024 Merritt Manufacturing, Inc. All Rights Reserved



**THE JEFFERSON**™  
Car seat specially designed for children with an  
omphalocele.



# FEATURES

## YOKE HARNESS



Easy to use yoke harness safely restrains the baby while avoiding contact with the omphalocele.

## SUPPORT CUSHIONS



Each Jefferson car seat comes with two support cushions, one for each buckle position, to optimally position smaller children.

## LOWER HARNESS ADJUSTERS



Lower harness adjusters ensure that the legs are properly secure and the lower harness is tightened while still avoiding contact with the omphalocele.

