

Dual Access Sling

- Provides support from under thighs, behind the mid-to-low back, and around the torso with a hook and loop belt enclosure
- For use with patients who are seated requiring transfer to tub or toileting
- Optional side retention straps
- Divided, padded leg supports that can be applied and removed while in a seated position

8B0X00




Product Codes and Size Chart:



Size	Product No.	A	C	D	E	F	G	H	I	J	K	L	M	N	Max SWL	Size Trim
XXS	8B0700	25" (64cm)	8½" (22cm)	9½" (24cm)	26¾" (68cm)	13" (33cm)	8¼" (21cm)	8" (20cm)	13" (33cm)	6¼" (16cm)	4¾" (12cm)	10½" (27cm)	18½" (47cm)	6" (15cm)	450lbs	Silver
XS	8B0600	28¾" (73cm)	8½" (22cm)	10¼" (26cm)	30" (76cm)	13¾" (35cm)	14" (36cm)	8" (20cm)	13¾" (35cm)	7" (18cm)	4¾" (12cm)	10½" (27cm)	18½" (47cm)	6" (15cm)	450lbs	Red
S	8B0500	31½" (80cm)	8½" (22cm)	11¾" (30cm)	34" (86cm)	17¼" (44cm)	14" (36cm)	10" (25cm)	17¼" (44cm)	8" (20cm)	4¾" (12cm)	10½" (27cm)	22" (56cm)	8" (20cm)	625lbs	Yellow
M	8B0400	37" (94cm)	8½" (22cm)	11¾" (30cm)	38" (97cm)	20½" (52cm)	14" (36cm)	10" (25cm)	20½" (51cm)	12" (30cm)	5½" (14cm)	10½" (27cm)	22" (56cm)	8" (20cm)	625lbs	Green
L	8B0300	42" (107cm)	10½" (27cm)	13¼" (34cm)	46" (117cm)	21½" (55cm)	14" (36cm)	10" (25cm)	21½" (55cm)	13¼" (34cm)	6" (15cm)	10½" (27cm)	22" (56cm)	8" (20cm)	625lbs	Blue
XL	8B0200	46" (117cm)	10½" (27cm)	16½" (42cm)	52¾" (134cm)	24" (61cm)	15¼" (39cm)	10" (25cm)	24" (61cm)	15¼" (39cm)	6" (15cm)	10½" (27cm)	22" (56cm)	8" (10cm)	1000lbs	Black

Specifications:



Applications:

- From bed to commode
- From wheelchair to toilet
- Transfer to tub