

**350 lb. Weight Capacity**  
 Effective January 26th, 2026

**Bill To:**

Date: \_\_\_\_\_  
 P.O.#: \_\_\_\_\_  
 Quote#: \_\_\_\_\_  
 Account#: \_\_\_\_\_  
 Ordered By: \_\_\_\_\_  
 Supplier Company: \_\_\_\_\_  
 ATP (REQUIRED): \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 City, State, Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Email: \_\_\_\_\_

**Ship To:**

Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City/State: \_\_\_\_\_  
 Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Tag For: \_\_\_\_\_



**AltimateMedical**  
 creating products · changing lives

**ActiveAid BCB Bath Commode Bench**
**ActiveAid BCB Bath Commode Bench Standard Configuration Includes:**

A modular 18"W x 18"D height adjustable frame.

**Frame- must select one**

☒ PH5000 18"Wx18"D Black Stainless Steel Frame **\$750**

**Legs- must select one**


☐ PH5011 Legs w/Crutch Tips (4) (height seat range from floor to bottom frame-16.75"-19.75") **\$200**  
☐ PH5012 Legs w/Suction Cups (4) (height seat range from floor to bottom frame-16.75"-19.75") **\$245**  
☐ PH5013 Suction Cups (2), Crutch Tips (2) **\$220**

**Seats- must select one**


☐ PH5901 18"Wx18"D Front Open Contoured Seat- Gel (2.25" gel) **\$660**  
☐ PH5902 18"Wx18"D Closed Contoured Seat- Foam (2.25" foam) **\$560**  
☐ PH5903 18"Wx18"D Closed Contoured Seat- Gel (2.5" gel) **\$760**  
☐ PH5904 18"Wx18"D Front Open Contoured Seat- Foam (2.5" foam) **\$460**  
☐ PH5905 18"Wx18"D Left Side Opening Contoured Seat- Foam (2.25" foam) **\$560**  
☐ PH5906 18"Wx18"D Right Side Opening Contoured Seat- Foam (2.25" foam) **\$560**  
☐ PH5911 18"Wx18"D Left Side Opening Contoured Seat- Gel (2.5" gel) **\$760**  
☐ PH5912 18"Wx18"D Right Side Opening Contoured Seat- Gel (2.5" gel) **\$760**  
☐ PH5907 18"Wx18"D Rear Opening Contoured Seat- Foam (2.25" foam) **\$560**  
☐ PH5908 18"Wx18"D Solid- Foam (2.25" foam) **\$560**




Front Opening



Closed



Left Side Opening



Right Side Opening



Rear Opening

**Back**

☐ PH5101 Solid Contoured Foam Back **N/C**  
☐ PH5102 Sling Back **\$245**

**Handles- back required**

☐ PH5110 No Push Handle **N/C**  
☐ Push Handle Fixed (pair) (standard) **N/C**  
☐ Fold-Down Push Handle (pair) **\$200**

**Armrests**

☐ PH5200 Flip-Up Armrest, Height Adjustable-Left **N/C**  
☐ PH5201 Flip-Up Armrest, Height Adjustable-Right **N/C**  
☐ PH5210 2-Point Locking Armrest-Left **\$150**  
☐ PH5211 2-Point Locking Armrest-Right **\$150**  
 (a 2-point locking armrest and transfer bench cannot be used on the same side of the chair at the same time)

**Belts**

☐ PCW60884 Waist Belt (66" circumference with 12" overlap) **\$80**  
☐ PCT60884 Trunk Belt (66" circumference with 12" overlap) **\$80**  
☐ PC5310 Bodypoint® Trunk Belt **\$224**  
☐ PC5312 Bodypoint® Trunk Belt & Waist Belt **\$277**  
☐ PC5314 Bodypoint® Belt w/Quick Release Knobs **\$260**

**Benches- 1" lateral adjustment, a 2-point locking armrest and transfer bench cannot be used on the same side of the chair at the same time**

☐ PH5803 Floor Supported Transfer Bench-9"W (Narrow w/Crutch Tip Legs) **\$379**  
☐ PH5804 Floor Supported Transfer Bench-9"W (Narrow w/Suction Cup Legs) **\$400**  
☐ PH5800 Floor Supported Transfer Bench-14"W (Wide w/Crutch Tip Legs) **\$495**  
☐ PH5802 Floor Supported Transfer Bench-14"W (Wide w/Suction Cup Legs) **\$516**  
☐ PH5801 Tub Supported Short Transfer Bench-9"W (Narrow w/Quad Tip Legs) **\$379**

**Pans- (may limit access over toilets)**

☐ PC5702 Pan, Cover, and Hanger 3" Drop **\$210**

MATRIX

**A. Overall Frame Width**

Overall frame width is measured from the leftmost point of the frame to the rightmost point of the frame.

**B. Overall Frame Depth**



Overall frame depth is measured from the most rearward point of the frame to the most forward point of the frame.

**C. Seat to Floor Height**

Seat to floor height is measured from the floor to the top front edge of the standard 2" seat cushion.

**D. Under Seat Clearance**

Under Seat Clearance refers to the vertical distance from the lowest point beneath the seat cushion to the floor. This measurement represents the space available from the bottom of the commode seat frame to the floor, ensuring sufficient clearance to accommodate a toilet.



- A. Floor to top of toilet bowl \_\_\_\_"
- B. Floor to top of toilet seat - down \_\_\_\_"
- C. Floor to top of tank \_\_\_\_"
- D. Wall to front of toilet \_\_\_\_"
- E. Front of tank to front of bowl \_\_\_\_"
- F. Wall to front of tank \_\_\_\_"
- G. Width of bowl \_\_\_\_"
- H. Width of tank \_\_\_\_"
- IL. Nearest obstruction - left side \_\_\_\_"
- IR. Nearest obstructions - right side \_\_\_\_"

Description	Specifications	
	Imperial	Metric
Weight limit	350 lbs	158.8 kg
Seat to floor height	19.5 - 22.5 in	49.5 - 57.1 cm
Overall frame width	20.25 in	51.4 cm
Overall width with narrow bench	29 - 30 in	73.7 - 76.2 cm
Overall width with wide bench	34 - 35 in	86.4 - 88.9 cm
Overall frame length (front to back)	18.25 in	46.4 cm
Overall length with backrest	21.25 in	54.0 cm
Overall length with sling back rigidizer	24 in	61.0 cm
Overall length with fixed handles	26.75 in	67.9 cm
Overall length with fold-down handles	23.75 - 27.5 in	60.3 - 69.9 cm
Seat dimensions	18 x 18 in	45.7 x 45.7 cm
Narrow bench dimensions	20 x 9 in	50.8 x 22.9 cm
Wide bench dimensions	20 x 14 in	50.8 x 35.6 cm
Distance between bench cushion and seat cushion	0.5 - 1.5 in	1.3 - 3.8 cm
Transfer bench to floor height	19.5 - 23.5 in	74.9 - 59.7 cm
Under seat height clearance	16.75 - 19.75 in	42.5 - 50.2 cm
Under seat width clearance	14.5 in	36.8 cm
Tub wall clearance height	16.75 - 20.75 in	42.5 - 52.7 cm
Tub support bench height adjustment	2 in	5.1 cm
Top of tub support bench to bottom of tub support feet	6.5 - 8.5 in	16.5 - 21.6 cm
Adjustable armrest height from top of seat	8 - 14.5 in	20.3 - 36.8 cm
Fixed 2-point armrest height from top of seat	9 in	22.9 cm
Width between arm rests	18.25 in	46.4 cm
Back height from top of seat	12.5 - 16.5 in	31.8 - 41.9 cm
Back angle	7°	7°
Weight of frame & seat	22 lbs	10.0 kg
Weight of backrest	8.75 lbs	4.0 kg
Weight of narrow bench	9 lbs	4.1 kg
Weight of wide bench	12.75 lbs	5.8 kg
Floor to Crossbar Height	7.5 - 10.5 in	19.1 - 26.7 cm
Outside of Crutch Tip to Outside of Crutch Tip Width	17.5 in	44.5 cm