

Hammock Sling PolySlip with Head Support

8D1X20

- Provides support around the hips, sacral area, and under the thighs up to the shoulders with optional head support
- Divided leg supports that can be applied or removed in seated or supine positions
- Multiple loop and strap options for seated and semi-reclined positions
- Legs strap application in various positions for multi-level support
- Built-in handles for ultimate caregiver movement control




Product Codes and Size Chart:



Polyester SlipFit (PolySlip)

Size	Part No.	A	B	C	W	LF	L	S	HS	Max SWL	Size Trim Color
XXS	8D1720	27½" (70cm)	4½" (11cm)	4" (10cm)	25" (64cm)	19" (48cm)	23½" (60cm)	15" (38cm)	6½" (17cm)	450lbs	White
XS	8D1620	33½" (85cm)	6" (15cm)	5" (13cm)	29½" (75cm)	21½" (55cm)	23½" (60cm)	15" (38cm)	8" (20cm)	450lbs	Silver
S	8D1520	37½" (95cm)	6" (15cm)	6" (15cm)	33½" (85cm)	27" (69cm)	23½" (60cm)	15" (38cm)	10½" (27cm)	625lbs	Red
M	8D1420	42" (107cm)	7" (18cm)	7" (18cm)	38" (97cm)	28¾" (73cm)	23½" (60cm)	15" (38cm)	13¾" (35cm)	625lbs	Yellow
L	8D1320	43" (109cm)	8" (20cm)	8" (20cm)	41" (104cm)	32" (81cm)	23½" (60cm)	15" (38cm)	13¾" (35cm)	625lbs	Green
XL	8D1220	46½" (118cm)	8" (20cm)	9" (23cm)	46½" (118cm)	35" (89cm)	23½" (60cm)	15" (38cm)	13¾" (35cm)	1000lbs	Blue

Specifications:



Applications:

