


SILVER CALCIUM ALGINATE DRESSINGS

Vitále® STRONG™ Silver Calcium Alginate Dressings are meticulously engineered to deliver superior wound care through a unique blend of calcium alginate, sodium alginate, and sodium silver phosphate. These dressings form a gel-like covering upon contact with wound exudate, maintaining a moist environment conducive to healing and preventing the dressing from adhering to the wound. The addition of silver provides broad-spectrum antibacterial protection, effectively targeting a wide range of pathogens including *Staphylococcus aureus*, *Escherichia coli*, *Enterococcus faecalis*, *Pseudomonas aeruginosa*, *Klebsiella pneumoniae*, and *Staphylococcus epidermidis*. Sterilized by gamma radiation, these dressings are designed for single use, ensuring a safe and hygienic application for moderate to heavily exuding wounds.



USAGE

Type: Primary

Stages: II-IV

Drainage/Exudate: Light to Heavy

Wound Types:

Pressure Sores

Venous Ulcers

Arterial Ulcers

Diabetic Ulcers

Grafts & Skin Donor Sites

Surgical/Post-Op

Traumatic Wounds

FEATURES & BENEFITS

- **Broad-Spectrum Antibacterial Protection:** The silver component provides effective antibacterial action against a wide range of pathogens, reducing the risk of infection.
- **Optimal Moist Healing Environment:** The calcium alginate forms a gel upon contact with exudate, maintaining a moist wound environment that supports natural healing processes.
- **High Absorbency:** Efficiently absorbs and retains exudate, managing wound moisture levels and reducing the risk of maceration.
- **Non-Adherent:** The gel-like covering prevents the dressing from sticking to the wound, allowing for painless removal and minimizing trauma to the wound bed.
- **Versatile Application:** Suitable for a variety of wounds, including leg ulcers, pressure ulcers, diabetic ulcers, trauma wounds, postoperative wounds, and graft and donor sites.

Vitále® STRONG™ Silver Calcium Alginate Dressings are trusted by healthcare professionals for their advanced wound management capabilities. The combination of high absorbency and broad-spectrum antibacterial protection ensures that wounds are kept in an optimal moist environment, promoting faster healing and reducing the risk of infection. These dressings are versatile, easy to apply, and suitable for a wide range of wound types, making them an indispensable tool in both clinical and home care settings. By choosing Vitále® STRONG™ Silver Calcium Alginate Dressings, healthcare providers can deliver superior wound care that enhances patient outcomes and satisfaction.

ORDERING INFORMATION

CATALOG NUMBER	DESCRIPTION	SIZE	UNITS/BOX	HCPC CODE
20550	Vitále® STRONG™ Silver Calcium Alginate Rope	3/4" x 12" (2 cm x 30 cm)	5	A6199
20552	Vitále® STRONG™ Silver Calcium Alginate Dressing	2" x 2" (5 cm x 5 cm)	10	A6196
20554	Vitále® STRONG™ Silver Calcium Alginate Dressing	4" x 4" (10 cm x 10 cm)	10	A6196
20555	Vitále® STRONG™ Silver Calcium Alginate Dressing	4" x 5" (10 cm x 12.5 cm)	10	A6197
20556	Vitále® STRONG™ Silver Calcium Alginate Dressing	6" x 6" (15 cm x 15 cm)	5	A6197
20558	Vitále® STRONG™ Silver Calcium Alginate Dressing	4" x 8" (10 cm x 20 cm)	5	A6197