

Protekt® Heavy Duty Bariatric Rollator

The **Protekt® Bariatric Rollator** provides stability and independence for larger individuals. The frame is constructed with extra strength reinforced steel and has large reinforced forks. The 8" wheels provide improved maneuverability and safety. Equipped with easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Easily folds for storage and transport with built-in fold strap. Supports patient weights up to 500 pounds.



Easy-to-operate
Ergonomic
Hand Brakes



8" Wheels



Water-Resistant
Nylon Pouch

FEATURES

- Durable reinforced steel frame
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- 8" black casters with soft grip tires, ideal for both indoor and outdoor use

COLOR

Black
Finish



SPECIFICATIONS

Overall Width:	29"
Width b/w Hand Grips:	25"
Seat Dimensions:	18"(W) X 13.5"(D)
Seat-to-Floor Height:	22"
Weight Capacity:	500 lbs.
Adjustable Handle Height:	34" - 38"
Wheel Size:	8"
Product Weight:	24 lbs.
Carton Dimensions:	36" x 23" x 10"
Shipping Weight:	28.5 lbs.
Warranty:	Limited Lifetime

Item #	Description	Color	HCPSC Code
RLBAST-BLK	Protekt Rollator, Bariatric Steel	Black	E0149 + E0156