


Pull your patients into the comfort zone



Introducing

Hydrofera Blue ComfortCel® Interface Antibacterial Foam Dressing



Find comfort with *blue.*
www.hydrofera.com

Hydrofera**BLUE**®

Hydrofera Blue ComfortCel® Interface



The newest addition to our portfolio of Advanced Wound Care Dressings:

- Pre-moistened and ready to use
- Effectively absorbs and transfers bacteria-laden exudate away from the wound bed
- Non-cytotoxic, does not inhibit growth factors or fibroblasts while providing broad-spectrum antibacterial protection

Soft PVA foam is designed for patient comfort during dressing changes:

- Minimizes adherence to wound bed
- Excellent wet tensile strength for easy removal from packed tunnels, undermining and deep wounds

Ordering Information

Product No	Product Size	Qty/Box
HBCC4550	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	5
HBCC6675	6" x 6" x 0.75" (15.2 cm x 15.2 cm x 1.9 cm)	5

Indications for use:

- Pressure Ulcers
- Diabetic Ulcers
- Venous Stasis Ulcers
- Arterial Ulcers
- Superficial Burns
- Donor Sites
- Post Surgical Incisions
- Traumatic Wounds
- Skin Tears
- Abrasions
- Lacerations

For product questions and samples, please call
Customer Service at 1-860-337-7730.

Hydrofera Dressings manufactured by:
Hydrofera, LLC, 340 Progress Drive
Manchester, CT 06042 USA
www.hydrofera.com

Hydrofera Blue and the Hydrofera Blue logo are trademarks of Hydrofera, LLC. ©2022
64799 RN L22
MADE IN USA

Hydrofera®