

PRODUCT SPECIFICATIONS

Core Products® Dual Comfort CorPak™ Hot or Cold Therapy Pack





Table of Contents:

Instructions	2
Features/Benefits	3
Size Small	4
Size Medium	6
Size Large	8
Size Cervical	10
Size Tri-Sectional	12
Size Wrist Wrap	14



Core Products
INTERNATIONAL, INC.

HEATING INSTRUCTIONS

For 3x5 inch pack, heat in microwave for 15 seconds on high, then at 5 second intervals until desired temperature is reached. For larger packs, heat in microwave for 30 seconds on high, then at 10 second intervals until desired temperature is reached. Knead after each heating. Overheating may cause injury or cause pack to rupture and leak.



COOLING INSTRUCTIONS

For convenience, store the pack in freezer. Place pack in the freezer for a minimum of 1 hour prior to use.

Core Products® Dual Comfort CorPak™

- Two surfaces for controlled therapy to ease pain
 - Strapping system provides compression and holds pack in place
 - Nontoxic, biodegradable fill maintains therapeutic temperature for 20 minutes
-
- Deux surfaces pour un traitement contrôlé visant à soulager la douleur
 - Système de bandage qui permet d'exercer une compression et de maintenir le sachet en place
 - Rembourrage non toxique et biodégradable qui maintient une température thérapeutique pendant 20 minutes
-
- Dos superficies para una terapia controlada para aliviar el dolor
 - El sistema con correas de sujeción brinda compresión y mantiene la almohadilla en su lugar
 - El relleno no tóxico, biodegradable, mantiene la temperatura terapéutica por 20 minutos

SMALL - 3" X 5"



Core Products® Dual Comfort CorPak™

- Product size: Small
- Material:
 - Food Grade Propylene Glycol USP
 - Bentonite clay
 - Water
 - Cover: Polyethylene / Polyester
- Product weight:
 - 0.25 lbs
- Color options:
 - Blue
- Latex:
 - This product is not made with natural rubber latex
- Product Origin:
 - Made in USA
- Brand:
 - Core Products®

MEDIUM - 6" X 10"

| 6" |



Core Products® Dual Comfort CorPak™

- Product size: Medium
- Material:
 - Food Grade Propylene Glycol USP
 - Bentonite clay
 - Water
 - Cover: Polyethylene / Polyester
- Product weight:
 - 1.25 lbs
- Color options:
 - Blue
- Latex:
 - This product is not made with natural rubber latex
- Product Origin:
 - Made in USA
- Brand:
 - Core Products®

LARGE - 10" X 13"

| 13" |

| 10" |



 **Core Products**
INTERNATIONAL, INC.
Dual Comfort CorPak
Hot & Cold Pack
www.coreproducts.com

Use hot or cold. For convenient cold therapy, store in freezer. For heat therapy, place pack in microwave and heat for 30 seconds. Continue heating at 10 second intervals until desired temperature is reached. Knead pack after each heating interval.

CAUTION: Microwave power levels vary.
Overheating may cause injury.

WARNING: Do not sit on or place excess weight on pack.

FROST FREE
Made in USA
21234JO

Core Products® Dual Comfort CorPak™

- Product size: Large
- Material:
 - Food Grade Propylene Glycol USP
 - Bentonite clay
 - Water
 - Cover: Polyethylene / Polyester
- Product weight:
 - 2.75 lbs
- Color options:
 - Blue
- Latex:
 - This product is not made with natural rubber latex
- Product Origin:
 - Made in USA
- Brand:
 - Core Products®

CERVICAL - 6" X 20"



Core Products® Dual Comfort CorPak™

- Product size: Cervical
- Material:
 - Food Grade Propylene Glycol USP
 - Bentonite clay
 - Water
 - Cover: Polyethylene / Polyester
- Product weight:
 - 1.75 lbs
- Color options:
 - Blue
- Latex:
 - This product is not made with natural rubber latex
- Product Origin:
 - Made in USA
- Brand:
 - Core Products®

TRI-SECTIONAL - 11" X 15"



Core Products® Dual Comfort CorPak™


- Product size: Tri-Sectional
- Material:
 - Food Grade Propylene Glycol USP
 - Bentonite clay
 - Water
 - Cover: Polyethylene / Polyester
- Product weight:
 - 3.75 lbs
- Color options:
 - Blue
- Latex:
 - This product is not made with natural rubber latex
- Product Origin:
 - Made in USA
- Brand:
 - Core Products®

DUAL COMFORT CORPAK™ - WRIST WRAP

| 14" |



|
8 1/2"
|




|
2"
|

Strap Total Length: 28.5"

Core Products® Dual Comfort CorPak™

- Product size: Wrist Wrap
- Material:
 - Food Grade Propylene Glycol USP
 - Bentonite clay
 - Water
 - Cover: Polyethylene / Polyester
- Product weight:
 - 2.19 lbs
- Color options:
 - Blue
- Latex:
 - This product is not made with natural rubber latex
- Product Origin:
 - Made in USA
- Brand:
 - Core Products®



Core Products International, Inc. | 808 Prospect Avenue | Osceola, WI 54020
800-365-3047 | www.coreproducts.com