



# Raze

MSRP – USA – USD  
Effective Nov 15, 2025

## Raised Toilet Seat / Transfer Tub Bench

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

Raze modular frame for over-toilet use / tub transfer bench when combined with the optional transfer surface.  
User weight capacity 352 lbs.

### 1. Frame Configurations

Z4441	Raze Frame with tall leg set – for over-toilet / shower stall use	\$995
Z4442	Raze Frame with short leg set – for tub use <b>MUST SELECT TRANSFER SURFACE BELOW</b>	\$1,025
<b>Transfer Surfaces (for tub and shower stall use) – MUST select a Transfer Surface for use in tub (Z4442).</b> Optional for shower stall use.		
ZRZTSN	Transfer Surface – Narrow (8.5" x 16"D) [Suitable for most tubs]	\$495
ZRZTSW	Transfer Surface – Wide (12.5" x 16"D)	\$595

### 2. Raze Travel Kit


ZRAZETK	Raze Travel Kit – Includes travel case, protective foam layer, seat and accessories pouches	\$395
---------	---	-------



Z4441



Z4441 with optional  
ZRB1815 and ZCHL



Z4442 with optional ZRB1815 and ZRZTSN

### 3. Optional Components

ZRAB1815	Raze Back - 18"W x 15"H w/ wrap-around removable arm supports, molded vinyl back upholstery (Recommended for tub use)	\$425
ZUZR1815	Back Upholstery Tension-Adjustable 18"W x 15"H w/ band clamps and upholstery support clips for wrap-around back style (must order with ZRAB1815)	\$215
ZCHL	Removable Cane Handle – Left (Required when Chest Belt is ordered with molded upholstery)	\$65
ZCHR	Removable Cane Handle – Right (Required when Chest Belt is ordered with molded upholstery)	\$65
<b>Additional Leg Sets (for use over toilet, in shower stall and in tub)</b>		
ZRZSL	Short Leg Set – for tub use (when Z4441 is ordered)	\$125
ZRZTL	Tall Leg Set – for toilet use / shower stall use with Transfer Surface (when Z4442 is ordered)	\$95

**4. Seats – Standard** – Please note that the Molded and VFI seats, formerly designated 16"D and 19"D, are now designated 17"D and 20"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18"W x 17"D – Front Bridged	\$425
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$425
ZMS20FRO	18"W x 20"D – Front Access	\$135	ZVFI17LTO	18"W x 17"D – Left Access	\$425
ZMS20BRI	18"W x 20"D – Front Bridged	\$135	ZVFI17RTO	18"W x 17"D – Right Access	\$425
<b>Specialty Seats</b>			ZVFI17RRO	18"W x 17"D – Rear Access	\$425
ZSOM	Shower-Only Seat – Molded 18"W x 16"-21"D	\$425	ZVFI20BRI	18"W x 20"D – Front Bridged	\$495
ZEJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$195	ZVFI20FRO	18"W x 20"D – Front Access	\$495
ZEJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	ZVFI20LTO	18"W x 20"D – Left Access	\$495
			ZVFI20RTO	18"W x 20"D – Right Access	\$495
			ZVFI20RRO	18"W x 20"D – Rear Access	\$495
<b>Seat Options</b>					
ZSS05	Spacers for Front Seat Brackets – 1/2" (pr) – for adding dump angle (based on 18"D seat)				\$85
ZSS10	Spacers for Front Seat Brackets – 1" (pr) – for adding dump angle (based on 18"D seat)				\$95

## 5. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. [See Custom Seat Assessment Form – Resources tab at razdesigninc.com](#)

ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below)	Select if both width and depth are ≤ 21"	\$795
ZCEJR	Custom E&J Replica Seat (enter specs below)	Select if both width and depth are ≤ 21"	\$995

### Custom Specifications (for above custom seats)

Seat Frame Width: 18"  
 Custom Seat Width: \_\_ 18" / 20"  
 Custom Seat Depth: \_\_ 17" / 18" / 19" / 20"  
 Access Opening: \_\_ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged / FAR – Front and Rear  
 Access Opening Width: \_\_ 2.5" / 3.0" / 3.5" / 4.0" / 4.5" / 5.0" (3.5"W is standard)  
 Aperture Width: \_\_ 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support, unless ordering Transfer Surface, which requires 18"W seat only.
- The Ischial Pelvic Alignment System (IPAS) allows 2" fore/aft adjustment of the seat (seat depth and aperture position) on all seats.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, etc.

## 6. Arm Supports and Accessories for Back Support ZRB1815

ZASWL	Arm Supports (pr) Wraparound – Long	Standard
ZASWS	Arm Supports (pr) Wraparound – Short	N/C
<b>Arm Support Pads</b>		
ZARMLL	Molded pad long 13" – Left (fits ZASWL)	Standard
ZARMLR	Molded pad long 13" – Right (fits ZASWL)	Standard
ZARMSL	Molded pad short <sup>1</sup> 10" – Left (fits ZASWS)	N/C
ZARMSR	Molded pad short <sup>1</sup> 10" – Right (fits ZASWS)	N/C
ZTAPLL	Tubular pad long 13" – Left (in lieu of molded pads, fits ZASWL)	N/C
ZTAPLR	Tubular pad long 13" – Right (in lieu of molded pads, fits ZASWL)	N/C
ZTAPSL	Tubular pad short 10" – Left (in lieu of molded pads, fits ZASWS)	N/C
ZTAPSR	Tubular pad short 10" – Right (in lieu of molded pads, fits ZASWS)	N/C
ZFAPLL	Molded Flat Armpad Long 14"L x 3.25"W – Left (fits ZASWL)	\$90
ZFAPLR	Molded Flat Armpad Long 14"L x 3.25"W – Right (fits ZASWL)	\$90
ZFAPSL	Molded Flat Armpad Short 10"L x 3"W – Left (fits ZASWS)	\$90
ZFAPSR	Molded Flat Armpad Short 10"L x 3"W – Right (fits ZASWS)	\$90
ZTRL	Molded Arm Trough 18.5"L x 4.5"W – Left (fits ZASWL)	\$110
ZTRR	Molded Arm Trough 18.5"L x 4.5"W – Right (fits ZASWL)	\$110
ZPAML	Pivoting Arm Mount <sup>2</sup> – Left (for ZFAPLL / ZTRL)	\$165
ZPAMR	Pivoting Arm Mount <sup>2</sup> – Right (for ZFAPLR / ZTRR)	\$165
ZAPSB18	Anterior Posture Support Bar (for ZARML only)	\$225

<sup>1</sup> Not compatible with ZFAP / ZTR / ZPAM

<sup>2</sup> Recommended when Lateral Thoracic Supports (TLTL/S) and Molded Arm Troughs (ZTR) are both ordered

## 7. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts (Requires Cane Handles)		
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$50	ZCBM	Chest Belt – Medium (back widths 16"-18", 1-piece w/ hook and loop closure)	\$45
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$50	ZCBL	Chest Belt – Large (back widths 18"-22", 1-piece w/ hook and loop closure)	\$45
ZPBICM	Infection Control Pelvic Belt <sup>1</sup> – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$75	ZCBICM	Infection Control Chest Belt <sup>1</sup> – Medium (back widths 16"-18", 1-piece w/ side-release buckle)	\$75
ZPBICL	Infection Control Pelvic Belt <sup>1</sup> – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$85	ZCBICL	Infection Control Chest Belt <sup>1</sup> – Large (back widths 18"-22", 1-piece w/ side-release buckle)	\$85
ZPBBPM	Pelvic Belt – Bodypoint Aeromesh – Med	\$215	ZCBBPM	Chest Belt – Bodypoint Aeromesh – Med	\$245
ZPBBPL	Pelvic Belt – Bodypoint Aeromesh – Large	\$215	ZCBBPL	Chest Belt – Bodypoint Aeromesh – Large	\$245

<sup>1</sup> Non-absorbent belts – recommended for infection-control environments

## 8. Positioning Pads for Back Support ZRB1815 – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

ZLTSWAP1L	Lateral Thoracic Support for Wraparound Arm w/ Small Pad – Left	\$325
ZLTSWAP1R	Lateral Thoracic Support for Wraparound Arm w/ Small Pad – Right	\$325
ZLTSWAP2L	Lateral Thoracic Support for Wraparound Arm w/ Large Pad – Left	\$325
ZLTSWAP2R	Lateral Thoracic Support for Wraparound Arm w/ Large Pad – Right	\$325

## 9. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125
ZSC17	Splash Curtain <sup>1</sup> for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$325
ZSC20	Splash Curtain <sup>1</sup> for 20"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$345
ZSCC	Custom Splash Curtain for Custom Seat – Raises seat 1/2"	\$345
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZCPAN	Commode Pan (with brackets)	N/C

<sup>1</sup> For Left, Right, Front, and Rear Open, call for order codes and pricing.

## 10. Credits and Omissions

Z915	Omit Seat	-\$45
------	-----------	-------

## 11. Notes:

SPECIFICATIONS	Raised Toilet Seat	Raised Toilet Seat w/ Optional Back Support	Transfer Tub Bench	Transfer Tub Bench w/ Optional Back Support
WEIGHT CAPACITY	352 lbs			
OVERALL WIDTH – w/ Narrow Transfer Surface (8.5") – w/ Wide Transfer Surface (12.5")	20"-23.75"	23.25"-23.75"	29.75"-31.25" 33.75"-35.25"	
OVERALL LENGTH (front to rear)	20"-24.25"			
BACK SUPPORT HEIGHT ABOVE SEAT	N/A	15"	N/A	15"
BACK SUPPORT WIDTH	N/A	18"	N/A	18"
ARM SUPPORT HEIGHT ABOVE SEAT	9" (w/ molded seat)			
SEAT FRAME WIDTH	19.5"			
SEAT DEPTH (middle of IPAS range)	17" or 20"			
CLEARANCE HEIGHT UNDER SEAT – 0.5" increments	16.5"-19.5" (w/ molded seat) * without brackets +0.75"		N/A	
UNDER-SEAT CLEARANCE WIDTH	11.25"-15.5"		N/A	
SEAT-TO-FLOOR HEIGHT	20.25"-23.25" 0.5" increments		N/A	
TRANSFER SURFACE-TO-FLOOR HEIGHT	N/A		21"-24" (to floor outside tub)	
TUB WALL CLEARANCE HEIGHT	N/A		16.75"-19.75" (to floor outside tub)	
PRODUCT WEIGHT	19 lbs	26 lbs	26.5 lbs	33.5 lbs
FRAME AND HARDWARE MATERIAL	Aluminum			
<sup>1</sup> Lifetime warranty against corrosion				