

Alpha Express, Flat Top, Treatment Table with Full Shelf

MODEL 3200-E

- Gray, powder-coated, all-steel frame and gray laminate shelf
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)



Black



Royal Blue

SPECIFICATIONS

3200-E27	72" L x 27" W x 31" H
3200-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Alpha Express, Flat Top, Treatment Table

MODEL 3100-E

- Gray, powder-coated, all-welded, steel frame
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)

SPECIFICATIONS

3100-E27	72" L x 27" W x 31" H
3100-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2



Black



Royal Blue

Alpha Express, Flat Top, Treatment Table with Drawers

MODEL 3130-E

- Gray, powder-coated, all-steel frame, gray laminate shelf, and two, easy-clean, plastic drawers with gray, drawer fronts
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)



Black



Royal Blue

SPECIFICATIONS

3130-E27	72" L x 27" W x 31" H
3130-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

Alpha Express, Flat Top, Treatment Table with Shelving

MODEL 3300-E

- Gray, powder-coated, all-steel frame and gray, laminate shelving
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)



Black



Royal Blue

SPECIFICATIONS

3300-E27	72" L x 27" W x 31" H
3300-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Alpha Express, Treatment Table

MODEL 3010-E

- Gray, powder-coated, all-steel frame, and adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)

SPECIFICATIONS

3010-E27	72" L x 27" W x 31" H
3010-E30	72" L x 30" W x 31" H



Black



Royal Blue

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

Alpha Express, Treatment Table with Full Shelf

MODEL 3020-E

- Gray, powder-coated, all-steel frame, gray laminate shelf and an adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)



Black



Royal Blue

SPECIFICATIONS

3020-E27 72" L x 27" W x 31" H

3020-E30 72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

Alpha Express, Treatment Table with Drawers

MODEL 3013-E

- Gray, powder-coated, all-steel frame, gray laminate shelf, two, easy-clean, plastic drawers with gray drawer fronts and an adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)



Black



Royal Blue

SPECIFICATIONS

3013-E27	72" L x 27" W x 31" H
3013-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Alpha Express, Treatment Table with Shelving

MODEL 3030-E

- Gray, powder-coated, all-steel frame, gray laminate shelving and an adjustable backrest
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Easy bolt-through assembly
- Paper dispenser included
- Adjustable leg levelers
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Available with 27" and 30" tops
- Heavyweight, knit-backed upholstery in popular Black or Royal Blue
- 450 lbs. load capacity under normal use (204 kg)

SPECIFICATIONS

3030-E27	72" L x 27" W x 31" H
3030-E30	72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2



Black



Royal Blue