

Easy Access Poly

- Provides good support due to the structured design
- Suitable for residents with poor or reduced tone
- Suitable for active residents
- Supportive trunk section wraps around the torso to provide full support to the body
- Full integral head support
- Enlarged commode aperture enable access for hygiene purposes
- Counterbalance straps provided allow for natural movement
- Wash temp 85°C / 185°F
- Tumble dry low heat (max 50°C)

M921532



Product Codes and Specifications:



100% Polyester		A	B	C	D	E
XS	M921532-XS	29¼ - 33" (74 - 84cm)	8 - 9½" (20 - 24cm)	23¾ - 26" (60 - 66cm)	29½" (75cm)	24" (61cm)
S	M921532-S	29¼ - 33" (74 - 84cm)	8 - 9½" (20 - 24cm)	26 - 29¼" (66 - 74cm)	29½" (75cm)	24" (61cm)
M	M921532-M	32¼ - 37" (82 - 94cm)	9½ - 11¾" (24 - 30cm)	29¼ - 31½" (74 - 80cm)	33½" (85cm)	26½" (67cm)
L	M921532-L	37 - 41¾" (94 - 106cm)	11 - 13½" (28 - 34cm)	31½ - 34¾" (80 - 88cm)	37½" (95cm)	28¾" (73cm)
XL	M921532-XL	41¾ - 47¼" (106 - 120cm)	11¾ - 14¼" (30 - 36cm)	34¾ - 39½" (88 - 100cm)	42½" (108cm)	30" (76cm)

