

SPAN PressureGuard® Air Therapy Surfaces Matrix

		Non-powered Air Therapy				Alternating/Lateral Rotation					Microclimate Management ["Low Air Loss" (LAL)]				Bariatric	
		Custom Care®	Span-Care® Convertible	Custom Care® Convertible	Custom Care® Convertible LAL	APM	APM2™	Safety Supreme	Span-Care® Convertible w/control unit	Custom Care® Convertible w/control unit	Protocol®	Easy Air®	Easy Air LR®	Custom Care® Convertible LAL w/control unit	Bariatric APM	Easy Air® XL
Therapy Modalities	Non-powered pressure redistribution	✓	✓	✓	✓											
	Convertible to powered via add-on control unit		✓	✓	✓											
	Alternating Pressure					✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Lateral Rotation						✓	✓		✓			✓	✓		
	Microclimate Management (LAL)									✓	✓	✓	✓	✓	✓	
	Powered Flotation						✓	✓	✓		✓	✓	✓	✓		✓
Frame Types	Versions to accommodate recessed deck frames	✓	✓	✓	✓				✓	✓			✓	✓		
Control Unit Features	Alarm: (V)isual, (A)udible					V	V, A	V, A	V, A	V, A	V, A	V	V	V	V, A	✓
	Timed Auto Firm						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Mattress Construction	Raised Perimeter							Standard			Available					
	Geo-Matt® top surface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Safety Edge™ bolster design	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Can maintain air in power outage; (A) = Automatically	✓ (A)	✓ (A)	✓ (A)	✓ (A)	✓	✓	✓	✓	✓ (A)	✓	✓ (A)	✓ (A)	✓ (A)	✓	✓ (A)
	Heel Slope™; (F) = Foam-only heel section	✓ (F)	✓ (F)	✓ (F)	✓ (F)	✓	✓	✓	✓	✓ (F)	✓	✓	✓	✓ (F)	✓	✓
	100% Latex free, 100% DHEP-free	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Cover	Non-stretch nylon with fluid-repellent coating					✓		✓							✓	
	Bi-directional stretch with Shear Transfer Zone™ design	✓	✓	✓			✓		✓	✓						
	Zippered air flow coverlet with Air Diffusion Matrix™ layer				✓						✓	✓	✓	✓		✓
General	Medicare Code	E0373	E0373	E0373	E0373	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277	E0277
	Therapeutic Weight Limit (lbs.)	500	500	500	500	350	350	350	400	350	400	500	500 (350 LR)	500 (350 LR)	750	1000
	Warranty: Mattress, non-pro-rated	5 yrs.	5 yrs.	5 yrs.	5 yrs.	18 mos.	18 mos.	18 mos.	5 yrs.	5 yrs.	18 mos.	2 yrs.	2 yrs.	5 yrs.	18 mos.	2 yrs.
	Warranty: Control unit, non-pro-rated					18 mos.	18 mos.	18 mos.	18 mos.	2 yrs.	18 mos.	2 yrs.	2 yrs.	2 yrs.	18 mos.	2 yrs.
	Meets Federal Fire Safety Regulations	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

Span-Care® Convertible



Custom Care® Convertible LAL



APM2™



Safety Supreme™



Protocol®



Custom Care® Convertible



Bariatric APM



Easy Air® LR



Easy Air® XL



Custom Care®



APM

