

THE ACTIVEAID

BCB

***BATH
COMMODE
BENCH***



ActiveAid



ActiveAid *BCB*

The ActiveAid BCB (Bath Commode Bench) is an all-new configurable bathing and toileting chair with seats that allow you to choose one of four directions for perineal access, along with multiple bench options to bridge tub walls.

The contoured seat features a cut-out that can be ordered with front, rear, or side opening and are available with molded foam or gel options. A solid seat is available for showering when bathing is the primary use. Independent, height-adjustable chair and bench legs help maintain a level seating surface in the shower or when bridging a tub ledge. With configurable seating, armrests, backs, patient positioning, and transfer bench options, the BCB adapts to a wide range of user needs and bathroom layouts. The ActiveAid BCB has a **weight capacity of 350lbs (158kg)**.



BCB FEATURES

CONFIGURABLE SEATING

The BCB is a modular system allowing users and clinicians flexibility to optimize accessibility while supporting transfer needs and daily hygiene activities. With molded foam, gel, and solid seating options, the BCB can be configured for commode use, bathing, or both—giving users options for toileting and hygiene routines.

SEATING OPTIONS

- Seat openings available for four directions of access (Front, rear, left and right open) over the toilet
- 10 seat options available, including Molded Gel and Molded Foam
- Solid seat option available to maximize seating surface when bathing or transfer flexibility is the priority



BENCHES

A selection transfer bench options expands the modular system to accommodate different bathroom environments and transfer requirements. Both narrow and wide bench pads are available with height-adjustable legs, along with a tub-supported bench option intended for spaces where permanent modifications are not possible.



SOLID CONTOURED FOAM



SLING BACK



SLING BACK

BACKS

The BCB's back supports provide posterior postural support.

- Solid Contoured Foam Back or Sling Back options available
- Sling Back features adjustable and tailored tension straps
- A horizontal rear crossbar is included to stabilize the back canes if ordered with a Sling Back



FLIP-UP



FLIP-UP



FLIP-UP AND
2-POINT LOCKING



2-POINT LOCKING

ARMRESTS

The BCB offers multiple armrest options that facilitate transfers, positioning, and user stability, allowing configuration based on functional needs and transfer technique.

- Flip-Up, Height-Adjustable Armrests or 2-Point Locking Armrests available on left or right side
- Can be configured with a combination of Flip-Up Armrest on one side and a 2-Point Locking Armrest on the other
- Flip-up Armrest brackets are height adjustable



PANS

The BCB accommodates both a commode pan system or over-toilet use. The BCB can be ordered with a removable, low-splash design pan when over-toilet use is not indicated.

- Pan, cover, and hanger system with a preset 3" (7 cm) drop
- Designed for easy placement and removal
- Pan features raised inner spine running front to back, designed to mitigate splashing



FOLD-DOWN HANDLE



FOLD-DOWN HANDLE



FIXED HANDLE

HANDLES

The BCB can be configured with optional handles attached to the back canes.

- May provide another point of contact for stability and balance
- Fixed or Fold-Down Push Handles available
- If no handles are selected, back canes are capped



BODYPOINT® TRUNK BELT



BODYPOINT® TRUNK BELT
AND WAIST BELT



TRUNK BELT



WAIST BELT

BELTS & POSITIONING SUPPORT

Positioning belts are available for additional waist and trunk support.

- Waist, trunk, and combination belt options available
- Bodypoint® belt options offered, including quick-release configurations
- Adjustable to support individual positioning needs



SUCTION CUP LEGS



CRUTCH TIP LEGS

BASE STABILITY

The BCB offers multiple leg tip configurations to accommodate different flooring conditions and environments.

- Crutch-tip legs, suction cup legs, or mixed configurations available
- Suction cups provide added grip on smooth, non-porous surfaces
- Chair and Bench legs are independently height-adjustable to accommodate some uneven surfaces and maintain a level seating surface



AltimateMedical
creating products · changing lives



ActiveAid



EasyStand



Therafin



Zing



VISIT ACTIVEAID.COM FOR:

- Comprehensive product specifications
- Assembly, fitting, and adjustment guides, videos and instructions for use
- Request a quote using online configuration
- Request a free product demonstration

ALTIMATEMEDICAL 262 W. 1ST ST., MORTON, MN • ALTIMATEMEDICAL.COM • 1-800-342-8968

AltimateMedical maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without prior notification. Check with AltimateMedical for latest information.
FORM BCBRO 011926 REV A Copyright © 2026 AltimateMedical, Inc. All rights reserved. U.S. & international patents pending. Printed in the U.S.A.