

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:



Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/5/2026


Rabbit Up



A (Total Height) _____
B (Armpit Height) _____
C (Top of Pelvic Height) _____
D (Buttock Height) _____
E (Chest Width) _____
F (Pelvic Width) _____
Age _____
Weight _____

Features, incl. in standard price


Stander Frame
Centre Bar
Brake on Castors
Angling Spindle




Mobile Stander Frame
Centre Bar
Driving Wheels
Angling Spindle
Closed Pushrim
0° Camber std. Angle
Spoke protectors
Central brake
Brake on Casters




Measurements




Size 1



Size 2



Size 3



Size 4

Product Specs as a standing frame

| | | | | |
|-------------------|-----------|-----------|-----------|-----------|
| Overall width (G) | 20¼" | 21¾" | 23¾" | 23¾" |
| Length (H) | 27¾" | 31¼" | 31¼" | 36½" |
| Total mass | 19.8 lbs. | 23.1 lbs. | 24.2 lbs. | 29.3 lbs. |

Product Specs as a mobile stander

| | w/26" wheels | w/32" wheels | w/32" wheels | w/36" wheels |
|-------------------|--------------|--------------|--------------|--------------|
| Overall width (G) | 25" | 26¾" | 28¾" | 28¼" |
| Length (H) | 24½" | 29½" | 29½" | 35½" |
| Total mass | 29.75 lbs. | 34.1 lbs. | 35.25 lbs. | 40.75 lbs. |

| | Size 1 | Size 2 | Size 3 | Size 4 |
|-------------------------------|---------------|---------------|---------------|---------------|
| Overall measurements | | | | |
| Width, inside (G1) | 12¾" | 14¾" | 16¾" | 16¾" |
| Overall height (I) * | 26 ¼" – 39 ¼" | 33 ¾" – 46 ¾" | 39" – 51 ½" | 47 ¼" – 60 ¼" |
| Height from foot plate (I1) * | 21 ¾" – 29 ¼" | 29 ½" – 42 ½" | 32 ½" – 46 ¾" | 37 ¾" – 56" |
| Center bar angle (a1) | 0 to +30° | 0 to +30° | 0 to +30° | 0 to +30° |
| Foot plate angle (a2) | +15° to -10° | +15° to -10° | +15° to -10° | +15° to -10° |
| Camber | 8° | 6° | 6° | 4° |
| Center bar length | 19¾" | 27½" | 31½" | 39¼" |

User Measurement Guide for Stander

| | | | | |
|-------------------------|---------------|----------------|----------------|----------------|
| Recommended user height | 30" - 39" | 39" - 54" | 43" - 58" | 48" - 71" |
| Max armpit height | 29" | 42" | 46" | 56" |
| User weight max | 88.1 lbs. | 110.2 lbs. | 132.2 lbs. | 154.3 lbs. |
| Recommended user age | 1 - 3 yr./old | 4 - 10 yr./old | 5 - 12 yr./old | 7 - 18 yr./old |

* Measured in upright position and with the foot plate in the lowest position. For smaller children it is possible to raise the foot plate.

Rabbit Stander Options

Make a selection of either the Prone Stander or the Mobile Prone Stander

Rabbit-Up Prone Stander



| | |
|---|--------|
| Rabbit-Up, Stander, Yellow, Size 1 (E0638) (8621220-00) | \$1648 |
| Rabbit-Up, Stander, Blue, Size 2 (E0638) (8622220-00) | \$1690 |
| Rabbit-Up, Stander, Grey, Size 3 (E0638) (8623220-00) | \$1731 |
| Rabbit-Up, Stander, Silver, Size 4 (E0638) (8624220-00) | \$1858 |

Rabbit Mobile Prone Stander



| | |
|---|--------|
| 22" Wheels, Yellow, Size 1 (E0642) (8621222-21) | \$3159 |
| 24" Wheels, Yellow, Size 1 (E0642) (8621223-21) | \$3159 |
| 26" Wheels, Yellow, Size 1 (E0642) (8621224-21) | \$3159 |
| 29" Wheels, Blue, Size 2 (E0642) (8622225-21) | \$3957 |
| 32" Wheels, Blue, Size 2 (E0642) (8622226-21) | \$3957 |
| 29" Wheels, Grey, Size 3 (E0642) (8623225-21) | \$4000 |
| 32" Wheels, Grey, Size 3 (E0642) (8623226-21) | \$4000 |
| 32" Wheels, Silver, Size 4 (E0642) (8624226-21) | \$4126 |
| 36" Wheels, Silver, Size 4 (E0642) (8624227-21) | \$4253 |

Chest Supports

Supports and Pad



Chest Support and Pad

| | |
|--|-------|
| Chest Support and Pad, W: 7.75-10.75", Size 1 (861704) | \$246 |
| Chest Support and Pad, W: 7.75-10.75", Size 2 and 3 (861702) | \$246 |
| Chest Support and Pad, W: 9.75-12.5", Size 4 (861703) | \$271 |



Filler Pad

Not required for Surrounding Chest Support

| | |
|----------------------------|------|
| Filler Pad, Small (861421) | \$98 |
| Filler Pad, Large (861422) | \$98 |



Extension for Center Bar

Available with Chest and Surrounding Chest Supports

| | |
|-----------------------------------|------|
| Extension for Center Bar (861480) | \$98 |
|-----------------------------------|------|

Surrounding Chest Support and Options



Surrounding Chest Support

| | |
|--|-------|
| Surrounding Chest Support, C:18.5-25.25", Small (861810) | \$589 |
| Surrounding Chest Support, C: 26.5-30.75", Medium (861820) | \$613 |
| Surrounding Chest Support, C: 34.25-44.75", Large (861830) | \$649 |



Cover for Surrounding Chest Support

| | |
|--|-------|
| Cover for Surrounding Chest Support, C:18.5-25.25", Small (869810-160) | \$139 |
| Cover for Surrounding Chest Support, C: 26.5-30.75", Medium (869820-260) | \$171 |
| Cover for Surrounding Chest Support, C: 34.25-44.75", Large (869830-360) | \$182 |



Extra Pad for Surrounding Support

| | |
|---|------|
| Extra Pad for Surrounding Support, H: 4.75" x W: 2.25" (869225) | \$75 |
|---|------|

Trays



Trays

| | |
|--|-------|
| Tray, Plastic, Grey (862515) | \$463 |
| Tray With Storage Box (862501) | \$475 |
| Plastic Tray with Bowl and Cover (862525) | \$672 |
| <i>Mounting Hardware required for Tray with Bowl</i> | |

Tray Accessories

| | |
|--|-------|
| Height Adjustment of Tray (862510) | \$271 |
| <i>Required with Tray with Storage Box</i> | |

Mounting Hardware

Required Selection for Tray with Bowl

| | |
|--|-------|
| Removeable Hardware Kit for Low Mount of Tray on Center Bar (862525-Kit) | \$320 |
|--|-------|

Abdominal Support



Abdominal Supports

| | |
|-----------------------------------|-------|
| Abdominal Support (Each) (861411) | \$193 |
|-----------------------------------|-------|

Head Supports

Head Support, Head Support Cover, Connecting System for Head Support, Headrest Mounting Hardware, and Hip Support w/ Pommel are required for a complete Head Support

Head Support



| | |
|---|-------|
| Anatomic Head Support, Small (5099001) | \$285 |
| Anatomic Head Support, Medium (5099002) | \$285 |
| Anatomic Head Support, Large (5099003) | \$285 |

Connecting System for Head Support




| | |
|---|-------|
| Connecting System for Head Support, Size 1 and 2 (861424-3) | \$310 |
| Connecting System for Head Support, Size 3 and 4 (861424-4) | \$310 |

Head Support Cover


| | |
|---|-----|
| Anatomic Head Support Cover, Small (9923-83) | N/C |
| Anatomic Head Support Cover, Medium (9945-83) | N/C |
| Anatomic Head Support Cover, Large (9967-83) | N/C |

Headrest Mounting Hardware



| | |
|--|-------|
| Headrest Mounting Hardware, Short (98619-25) | \$250 |
| Headrest Mounting Hardware, Long (98619-26) | \$256 |

Pelvic and Hip Control Options



Hip Support w/ Pommel


Front mount to center bar. Required with Head Support.

| | |
|--|-------|
| Hip Support w/ Pommel, Size 1 (861501) | \$607 |
| Hip Support w/ Pommel, Size 2 and 3 (861502) | \$607 |
| Hip Support w/ Pommel, Size 4 (861503) | \$672 |



Hip Support, With Rotation

| | |
|---|-------|
| Hip Support, With Rotation, Size 2 and 3 (861612) | \$426 |
| Hip Support, With Rotation, Size 4 (861613) | \$463 |



Hip Support, Front Mount

| | |
|---|-------|
| Hip Support, Front Mount, Size 1 (861601) | \$283 |
| Hip Support, Front Mount, Size 2 and 3 (861602) | \$319 |
| Hip Support, Front Mount, Size 4 (861603) | \$352 |



Hip Support, Surrounding

Front mount to center bar

| | |
|---|-------|
| Hip Support, Surrounding, Size 1 (861610) | \$667 |
| Hip Support, Surrounding, Size 2 and 3 (861620) | \$667 |
| Hip Support, Surrounding, Size 4 (861630) | \$679 |



Pommel With Back Support

| | |
|-----------------------------------|-------|
| Pommel With Back Support (861401) | \$524 |
|-----------------------------------|-------|

Knee Supports



Knee Support

| | |
|---|-------|
| Knee Supports, Size 1 (Pair) (861301) | \$336 |
| Knee Supports, Size 2 and 3 (Pair) (861302) | \$426 |
| Knee Supports, Size 4 (Pair) (861303) | \$426 |




Hinged Knee Supports

| | |
|--|-------|
| Hinged Knee Supports, Size 1 (Pair) (861311) | \$524 |
| Hinged Knee Supports, Size 2 and 3 (Pair) (861312) | \$524 |
| Hinged Knee Supports, Size 4 (Pair) (861313) | \$524 |

Footplate and Sandals Options

Flip-up Footplate required for all Sandals, except the Individual Footrest w/ Hinged Sandals

Footplates and Options




Flip-Up Footplate

| | |
|---|-------|
| Flip-Up Footplate, Size 1 and 2 (861211) | \$265 |
| Flip-Up Footplate, Size 2, 3 and 4 (861212) | \$265 |



Formed Footplates

| | |
|--|-------|
| Formed Footplates, Size 0 (Pair) (86149) | \$198 |
| Formed Footplates, Size 1 and 2 (Pair) (86147) | \$198 |
| Formed Footplates, Size 3 (Pair) (86148) | \$198 |



Footplate Straps

| | |
|--|------|
| Foot Straps, Plastic Lock (Pair) (85096) (K0040) | \$42 |
| Footplate Straps, Velcro (Pair) (85102) | \$59 |



Individual Footrest w/ Hinged Sandals

| | |
|---|-------|
| Individual Footrest w/ Hinged Sandals, Size 1 (Pair) (861201) | \$475 |
| Individual Footrest w/ Hinged Sandals, Size 2 and 3 (Pair) (861202) | \$475 |

Sandals and Options



Sandals

| | |
|---|-------|
| Sandals for Rabbit Up/ Toucan, Small (869490-RT) | \$233 |
| Sandals for Rabbit Up/ Toucan, Medium (869491-RT) | \$244 |
| Sandals for Rabbit Up/ Toucan, Large (869492-RT) | \$251 |



Sandals

| | |
|--|-------|
| Hinged Sandals w/ Straps, Size 1 (Pair) (861221) | \$162 |
| Hinged Sandals w/ Straps, Size 2 and 3 (Pair) (861222) | \$162 |
| Hinged Sandals w/ Straps, Size 4 (Pair) (861223) | \$283 |



Heel Stops

| | |
|-----------------------------------|------|
| Heel Stops, Small (Pair) (861251) | \$61 |
| Heel Stops, Large (Pair) (861252) | \$71 |

Conversion Kits

Please quote two for a complete set of wheels. Two wheels are required for a full conversion kit. Conversion kit is used to convert Prone Stander into Mobile Prone Stander.

Driving Wheel Conversion Kits



| | |
|---|--------|
| Driving Wheels 22", Closed Pushrims Size 1 (Each) (134209) | \$694 |
| Driving Wheels 24", Closed Pushrims Size 1 (Each) (134210) | \$694 |
| Driving Wheels 26", Closed Pushrims Size 1, 2 and 3 (Each) (134211) | \$694 |
| Driving Wheels 29", Closed Pushrims Size 2 and 3 (Each) (134212) | \$1071 |
| Driving Wheels 32", Closed Pushrims Size 2, 3 and 4 (Each) (134213) | \$1071 |
| Driving Wheels 36", Closed Pushrims Size 4 (Each) (134214) | \$1071 |



Central Brake System

| | |
|---|-------|
| Central Brake System, Size 1 (134215) | \$349 |
| Central Brake System, Size 2 and 3 (134216) | \$349 |
| Central Brake System, Size 4 (134218) | \$349 |



Camber Kit Rabbit Up

| | |
|--------------------------------------|------|
| Camber Kit Rabbit Up 0° (134200) | \$63 |
| Camber Kit Rabbit Up +3/-3° (134201) | \$63 |
| Camber Kit Rabbit Up +6/-6° (134202) | \$63 |