



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E

- Child-friendly design to encourage therapy compliance
- Convenient Doggy House Backpack case comes inside each Beagle Neb box
- Fun and entertaining doggy stickers included
- Complete with disposable neb kit, mouthpiece, pediatric aerosol mask, replacement filter and instructions
- Powerful piston pump is efficient, dependable and quiet
- Retail packaged

PRODUCTS

Description

18090-BE	Pediatric Beagle Nebulizers,1/ea, 6/cv
18091-BE	Ped.BeagleNebReuse&DispNebKit, 1/ea,6/cv

SPECIFICATIONS

Description

Electrical Requirements	120 VAC 60Hz
Liter Flow	8 lpm
Maximum Pressure	42 PSI
Method of Nebulization	Piston Compressor
Operating Pressure	13 PSI
Particle Size	0.5 µm to 10 µm
Unit Dimensions	5.6" (W) x 8.6" (H) x 8.6" (D)
Warranty	Limited 5 Year
Weight	3.1 lbs

ACCESSORIES

Description

18090F	NebFilter 10/bag, BAG
18031	Nebulizer Carry Bag, 1/ea, EA
18254	Reusable Neb Kit 10/case, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS