

Company: Ship To:

Account Number: Date: Address:










☐ Quote: ☐ PO Number: City: State:

ATP: Zip: Phone:

Marked For: Email/Fax:

Associated Quote/Order#: Effective: 1/5/2026

Measurements and Fitting

		Size 0	Size 1	Size 1 AT	Size 2	Size 2 AT	Size 3	Size 3 AT
	Width between handles	100-250 mm (3.9" - 9.8")	190-350 mm (7.5" - 13.8")	190-350 mm (7.5" - 13.8")	240-420 mm (9.4" - 16.5")	240-420 mm (9.4" - 16.5")	330-480 mm (13" - 18.9")	330-480 mm (13" - 18.9")
	Height from floor to handles	300-610 mm (11.8" - 24")	420-710 mm (16.5" - 28") 380-700 mm* (15" - 27.6")	460-750 mm (18.1" - 29.5")	560-850 mm (22" - 33.5")	600-890 mm (23.6" - 35")	740-1025 mm (29.1" - 40.4")	780-1065 mm (30.7" - 41.9")
	Max. user height	1050 mm (41.3")	1150 mm (45.3")	1190 mm (46.9")	1500 mm (59")	1540 mm (60.6")	1800 mm (70.9")	1840 mm (72.4")
	Max. user mass	30 kg (66 lbs)	30 kg (66 lbs)	30 kg (66 lbs)	45 kg (99 lbs)	45 kg (99 lbs)	80 kg (176 lbs)	80 kg (176 lbs)
	Length	640 mm (25.2")	720 mm (28.3")	780 mm (30.7")	840 mm (33.1")	900 mm (35.4")	870 mm (34.3")	930 mm (36.6")
	Overall width	580 mm (22.8")	630 mm (24.8")	660 mm (26")	690 mm (27.2")	700 mm (27.6")	690 mm (27.2")	720 mm (28.3")
	Folded length	650 mm (25.6")	720 mm (28.3")	750 mm (29.5")	860 mm (33.9")	890 mm (35")	1000 mm (39.4")	1030 mm (40.6")
	Folded width	580 mm (22.8")	620 mm (24.4")	660 mm (26")	690 mm (27.2")	700 mm (27.6")	690 mm (27.2")	720 mm (28.3")
	Folded height	300 mm (11.8")	300 mm (11.8")	345 mm (13.6")	270 mm (10.6")	310 mm (12.2")	290 mm (11.4")	330 mm (13")

*Mounted with the Starter kit accessory. Further measurements are available on our website.

Colors



Crocodile



Crocodile All-Terrain

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1197.5 Pricing and specifications subject to change without notice.

Crocodile Frame

Make a selection of either the standard Crocodile or the Crocodile AT

Crocodile Base

<input type="checkbox"/> Crocodile, White, Size 0 (86800)	\$733	
<input type="checkbox"/> Crocodile, Yellow, Size 1 (86801)		\$808
<input type="checkbox"/> Crocodile, Orange, Size 2 (86802)		\$808
<input type="checkbox"/> Crocodile, Metallic Black, Size 2 (86852)		\$808
<input type="checkbox"/> Crocodile, Anthracite Grey, Size 3 (86803)		\$1333

Crocodile AT Base

<input type="checkbox"/> Crocodile AT, Red, Size 1 (86801AT)	\$975	
<input type="checkbox"/> Crocodile AT, Green, Size 2 (86802AT)		\$1025
<input type="checkbox"/> Crocodile AT, Grey, Size 3 (86803AT)		\$1525

Anti-Tip Device



Swing-Away Anti-Tips for Crocodile

- ☐ Swing-Away Anti-Tips, for Croc Size 0-2 (86815) \$175
- ☐ Swing-Away Anti-Tips, for Croc Size 3 (86815-3S) \$185



Fixed Anti-Tips for Crocodile

- ☐ Anti-Tip Fixed, Pair, for Croc Size 0-2 (86815-F) \$85
- ☐ Anti-Tip Device Fixed, for Croc Size 3 (86815-3F) \$89

Swing-Away Anti-Tips for Crocodile AT

- ☐ Swing-Away Anti-Tips for Croc AT Size 1-2 (868151) \$225
- ☐ Swing-Away Anti-Tips for Croc AT Size 3 (868153) \$245

Brake



Hand Brake Kits

Eliminates the function of anti-reverse wheel stop

- ☐ Hand Brake Kit, for Size 1-2 (86834) \$270
- ☐ Hand Brake Kit, for Size 3 (86873) \$270



Deceleration Hand Brake Kit

Not available with Croc AT. Eliminates the function of anti-reverse wheel stop.

- ☐ Deceleration Hand Brake Kit, Sizes 1-2 (86834-DEC) \$270
- ☐ Deceleration Hand Brake Kit, Size 3 (86873-DEC) \$270

Drum Brake Kit

Allows for use of the anti-reverse wheel stop

- ☐ Drum Brake Kit, Size 1 (86834-1DB) \$995
- ☐ Drum Brake Kit, Size 2 (86834-2DB) \$995
- ☐ Drum Brake Kit, Size 3 (86873-3DB) \$995

Decelerating Hand Brake for Croc AT

- ☐ Deceleration Hand Brake for Croc AT, Size 1 (868341) \$295
- ☐ Deceleration Hand Brake for Croc AT, Size 2 (868342) \$295
- ☐ Deceleration Hand Brake for Croc AT, Size 3 (868343) \$295

Holder for Hand Brake

Required with use of Hand Brake and Forearm Support. Allows for mounting of Hand Brake to Forearm Support.

- ☐ Small, for Size 1 (86874-1) \$367
- ☐ Medium, for Size 2 (86874-2) \$414
- ☐ Large, for Size 3 (86874-3) \$445

Handles and Grip Bars

Grip Bar

Grip Bar is available for Crocodile sizes 0, 1, and 2



Standard Grip Bar

- ☐ Grip Bar, for Size 0 (86879-0) _____ \$141
- ☐ Grip Bar, for Size 1 (86879-1) _____ \$163
- ☐ Grip Bar, for Size 2 (86879-2) _____ \$196



Angle Adjustable Grip Bar

Angle Adjustable Grip Bar requires Surrounding Support or Back and Lateral Supports

- ☐ Angle Adjustable Grip Bar, Size 0 (86878-0) _____ \$235
- ☐ Angle Adjustable Grip Bar, Size 1 (86878-1) _____ \$275
- ☐ Angle Adjustable Grip Bar, Size 2 (86878-2) _____ \$305

Attendant Handle



Attendant Guide Bar

Available only for Crocodile sizes 0, 1 and 2

- ☐ Attendant Handle (86866-1) _____ \$295

Arm and Forearm Options

Arm Supports



Forearm Support with Handgrip

- ☐ Forearm Support Pair - for Sizes 0-1 (86858-1) _____ \$495
- ☐ Forearm Support Pair - for Size 2 (86859-1) _____ \$495
- ☐ Forearm Support Pair - for Size 3 (86860-1) _____ \$505



Handle Knobs

- ☐ Handle Knobs, Pair (86861) _____ \$56

Padded Forearm Positioner

For use with Forearm Support with Handgrip

- ☐ Padded Forearm Positioner Strap (869370) _____ \$75

Quantity:



Starter Kit for Crocodile

Available for Crocodile size 1



- ☐ Crocodile Starter Kit, for Croc size 1 (86801-0) _____ \$495

Back and Pelvic and Hip Supports



Surrounding Supports

- ☐ Surrounding Support for Size 0 (868610-0S) \$595
- ☐ Surrounding Support for Size 1 (868610-1S) \$595
- ☐ Surrounding Support for Size 2 (868610-2M) \$650



Back Support

- ☐ Back Support, Adjustable Padded for Size 3 (86886) \$552



Pelvic Support With Laterals

- ☐ Pelvic Support With Laterals, for Sizes 0-2 (86808) \$325
- ☐ Pelvic Support With Laterals, For Sizes 1-2 (86809) \$325

Abdominal Support Belt

For size 3 Back Support only

- ☐ Abdominal Support Belt, for Size 3 (85490-2) \$134



Lateral Strap

Used with Pelvic Support with Laterals

- ☐ Lateral Strap (SN4016) \$65



Mounting Hardware

Required for use of Pelvic and Back Support on size 3

- ☐ Cross Bar for Mounting Back/Pelvic Supports, for Size 3 (86880) \$107



Hip Support Pair

- ☐ Hip Supports, Pair (86819) \$335



Pad for Back Support

- ☐ Pad for Back Support (86816) \$195

Seating Options



Solid Flip Down Seat

- ☐ Size 0 (86820-01) \$186
- ☐ Size 1 (86820) \$186
- Not compatible with Crocodile starter kit*
- ☐ Size 2 (86822) \$186
- ☐ Size 3 (86863) \$189



Seat Sling

Not Available with Size 3

- ☐ Seat Sling, for Size 1 (86824) \$185
- ☐ Seat Sling, for Size 2 (86825) \$185



Lightweight Sling Seat

Requires Pelvic Support with Laterals and Lateral Strap or Surrounding Supports

- ☐ Lightweight Sling Seat, for Size 0-1 (86824-1) \$147
- ☐ Lightweight Sling Seat, for Size 2 (86825-1) \$147

Accessories



- ☐ Storage Bag for Sizes 1-2 (86839) \$78
- ☐ Storage Bag, for Size 3 (86843) \$230