



# Raz-AP

## Attendant Propel

MSRP – USA – USD  
Effective October 01, 2025

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

**Raz-AP:** Attendant Propel Mobile Shower Commode Chair with height-adjustable stainless-steel frame, Fixed Back, flip-up padded arm supports, foot supports, four dual-locking 5" (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. **User weight capacity 407 lbs.**

### 1. Standard Raz-AP Frame (18"W) – If you require a special width, depth and / or height see [Sec. 2]

Z10018	Raz-AP w/ 18"W base frame – Standard base frame width, will fit over most toilets	\$1,695
--------	---	---------

### 2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.7]

Z10016	Raz-AP w/ 16"W frame width <sup>1</sup> (base chair price – custom upcharge below will be added)	\$1,695
Z1SW16	Upcharge for custom 16"W frame – will not fit over most toilets	\$495
Z10020	Raz-AP w/ 20"W frame width <sup>2</sup> (base chair price – custom upcharge below will be added)	\$1,695
Z1SW20	Upcharge for custom 20"W frame	\$495
Z10022	Raz-AP w/ 22"W frame width (base chair price – custom upcharge below will be added)	\$1,695
Z1SW22	Upcharge for custom 22"W frame	\$595
<b>Custom Frame Configurations</b>		
ZLNGF3	Upcharge for custom 3" Longer Frame (3" longer at the front)	\$395
Z1LOWF4	Upcharge for custom Low Frame <sup>3</sup> (16" - 18" STF, 4" lower than standard)	\$395
Z1APRF3	Upcharge for Recliner Frame (3" longer at the rear for added stability)	\$395

<sup>1</sup> 18"W Visco-Foam Interface seats can be modified to fit.

<sup>2</sup> 18"W seats can be selected but must also select Seat Support Kit ZSSKAS [Sec.3]

<sup>3</sup> Forward Offset Foot Support Receiver ZFOR [Sec. 11] is required when selecting seats that are 20" or deeper.

**3. Seats – Standard** – Please note that the Molded and VFI seats, formerly designated 16"D and 19"D, are now designated 17"D and 20"D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Molded Seats (thickness ~2.5"-3")			Visco-Foam Interface Seats (thickness ~3.5")		
ZMS17FRO	18"W x 17"D – Front Access	Standard	ZVFI17BRI	18"W x 17"D – Front Bridged	\$425
ZMS17BRI	18"W x 17"D – Front Bridged	N/C	ZVFI17FRO	18"W x 17"D – Front Access	\$425
ZMS20FRO	18"W x 20"D – Front Access	\$135	ZVFI17LTO	18"W x 17"D – Left Access	\$425
ZMS20BRI	18"W x 20"D – Front Bridged	\$135	ZVFI17RTO	18"W x 17"D – Right Access	\$425
<b>Specialty Seats</b>			ZVFI17RRO	18"W x 17"D – Rear Access	\$425
ZSOM	Shower-Only Seat – Molded 18"W x 16"-21"D	\$425	ZVFI20BRI	18"W x 20"D – Front Bridged	\$495
ZEJRM	E&J Replica Seat – Molded – 18"W x 17"D	\$195	ZVFI20FRO	18"W x 20"D – Front Access	\$495
ZEJRV	E&J Replica Seat – Visco Foam – 18"W x 17"D	\$895	ZVFI20LTO	18"W x 20"D – Left Access	\$495
			ZVFI20RTO	18"W x 20"D – Right Access	\$495
			ZVFI20RRO	18"W x 20"D – Rear Access	\$495
<b>Seat Options</b>					
ZSSKAS	Seat Support Kit for Raz-AP – These are used to support the sides of seats that are narrower than seat frames on which they are mounted. A typical example would be when a 20"W chair is needed but an 18"W standard seat (Molded, VFI, EJR, or Shower Only Seat) is selected.				\$145
ZSS05	Spacers for Front Seat Brackets – 1/2" (pr) – for adding dump angle (based on 18"D seat)				\$85
ZSS10	Spacers for Front Seat Brackets – 1" (pr) – for adding dump angle (based on 18"D seat)				\$95

#### 4. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. [See Custom Seat Assessment Form](#) – Resources tab at razdesigninc.com

ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below) <b>Select if both width and depth are ≤ 21"</b>	\$795
ZCEJR	Custom E&J Replica Seat (enter specs below) <b>Select if both width and depth are ≤ 21"</b>	\$995
ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) <b>Select if width and/or depth is ≥ 22"</b>	\$1,095

#### Custom Specifications (for above custom seats)

Base Frame Width: \_\_\_ 16" / 18" / 20" / 22"  
 Custom Seat Width: \_\_\_ 16" / 18" / 20" / 22" / 24" (cannot be more than 2" wider than frame width)  
 Custom Seat Depth: \_\_\_ 15" / 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires Long Frame (Z1LNGF3)]  
 Access Opening: \_\_\_ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged  
 Access Opening Width: \_\_\_ 2.5" / 3.0" / **3.5"** / 4.0" / 4.5" / 5.0" (**3.5"W** is standard)  
 Aperture Width: \_\_\_ 3.0"-8.0" (6.5" recommended) – enter "NA" if aperture is not required

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seats.
- Custom Seat width should not be narrower than frame width.
- Custom Seat can be up to 2" wider than seat frame width to provide additional thigh support without ordering wider chair frame.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 1" fore / aft adjustment from specified custom seat depth.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.

#### 5. Fixed Back Frames (18"W only)

This is called a "Fixed" Back because the arm supports are not adjustable in height. Includes molded vinyl back support, fixed-height flip-up padded arm supports, back angle adjustment (90°/97°) with cane-style posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

ZFB1818	Fixed Back – 18"W x 18"H (Arm Support Lock Kit (ZALK) [Sec.9] is not included.)	Standard
ZFBU1818	Tension-Adjustable Back Upholstery for Fixed Back	\$185

## 6. Adjusta-Back Frames – Back frame widths will be the same as the seat frame widths unless otherwise selected in Sec 8

Includes height- and angle-adjustable arms with Armrest Lock Kit (ZALK), tension-adjustable back upholstery, flip-up padded arm supports, back angle adjustment (90°/97°). Back heights of 15"/17"/20" include cane-style posts, on which back upholstery can be positioned 2.5" higher than nominal height indicated below. e.g. 20"H can be moved up to 22.5"H. Adjusta-Back – 24"H includes straight posts, on which the upholstery cannot be raised.

ZAB __ 15*	Adjusta-Back <sup>1</sup> – 15"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	\$225
ZAB __ 17*	Adjusta-Back <sup>1</sup> – 17"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	\$225
ZAB __ 20*	Adjusta-Back <sup>1</sup> – 20"H	* insert back frame width to complete the code (14" / 16" / 18" / 20" / 22" / 24")	\$295
ZAB __ 24*	Adjusta-Back <sup>1</sup> – 24"H	* insert back frame width to complete the code (16" / 18" / 20" / 22" / 24")	\$345

**Reclining Backs – 30° (must select Reclining Frame option Z1APRF3 in Sec 2)** Reclining back frame width must be the same as the seat frame width. 2-Point Arms (Z2PLAS / Z2PSAS) are included with Reclining Backs.

ZRBA __ 20*	Reclining Adjusta-Back with 20"H back (anti-tippers included)	* insert back frame width to complete the code (16" / 18" / 20" / 22")	\$1,095
ZRBA __ 24*	Reclining Adjusta-Back with 24"H back (anti-tippers included)	* insert back frame width to complete the code (16" / 18" / 20" / 22")	\$1,195
ZDTL	Dual Triggers – both on left side		N/C
ZDTR	Dual Triggers – both on right side		N/C

<sup>1</sup> Dual Tilt Triggers option positions both triggers on one side, allowing for single-handed activation.

## 7. Custom Back Frame Width Narrower Than Chair Frame Width – Must choose back frame [Sec.7] (configurations below are not available for recliner backs)

**NOTE: For back frame width variances on recliner frames, please request custom quotation.**

ZNB__ __ *	Custom – Back Frame 2" Narrower than chair frame	\$345	Chair Frame Width	16"	18"	20"	22"	
			Back Frame Width	14" <sup>1</sup>	16" <sup>2</sup>	18" <sup>3</sup>	20" <sup>4</sup>	

\* Enter chair and back frame widths to complete code.  
eg. ZNB2018 specifies a 20"W chair frame with an 18"W back frame

<sup>1</sup> Requires Custom Seat

<sup>2</sup> Compatible with all 18"W seats but with reduced seat depth range. Also compatible with Custom or modified Visco-Foam Interface seats.

<sup>3</sup> Custom Seat recommended. Also compatible with all 18"W seats with Seat Support Kit Z4SSKAS.

<sup>4</sup> Custom Seat recommended.

## 8. Back Support Options – These back supports replace the molded vinyl upholstery on a Fixed Back or tension-adjustable upholstery on an Adjusta-Back. Symphony and Harmony Backs are recommended for infection-control environments.

ZSB14	14"W Symphony Back <sup>1</sup> – recessed solid back support	\$295
ZSB16	16"W Symphony Back <sup>1</sup> – recessed solid back support	\$295
ZSB18	18"W Symphony Back <sup>1</sup> – recessed solid back support	\$295
ZSB20	20"W Symphony Back <sup>1</sup> – recessed solid back support	\$395
ZSB22	22"W Symphony Back <sup>1</sup> – recessed solid back support	\$465
ZHB	Harmony Back <sup>1,2</sup> (contoured solid back – 18"W x 13"H) Fits 18"W back frames only Fits 17" / 20" back heights only	\$295

<sup>1</sup> Does not fit 15"H backs

<sup>2</sup> Not recommended for 24"H backs or reclining backs.

## 9. Arm Supports and Accessories

ZARMTL	Long Arm Support Tubes (pr)	Standard
ZARMTS	Short Arm Support Tubes (pr) <sup>1</sup>	N/C
ZARMLL	Long Arm Support Pad 13"L – Left	Standard
ZARMLR	Long Arm Support Pad 13"L – Right	Standard
ZARMSL	Short Arm Support Pad 10"L – Left	N/C
ZARMSR	Short Arm Support Pad 10"L – Right	N/C
ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	\$75
Z2PLAS	2-Point Arm Supports <sup>2</sup> (pr) with 13"L pads (std with reclining backs)	\$345
Z2PSAS	2-Point Arm Supports <sup>2</sup> (pr) with 10"L pads (std with reclining backs)	\$345
ZFAPLL	Long Molded Flat Armpad <sup>3</sup> – Left (14"L x 3.25"W – fits ZARM13)	\$90
ZFAPLR	Long Molded Flat Armpad <sup>3</sup> – Right (14"L x 3.25"W – fits ZARM13)	\$90
ZFAPSL	Short Molded Flat Armpad <sup>3</sup> – Left (10"L x 3"W – fits ZARM10)	\$90
ZFAPSR	Short Molded Flat Armpad <sup>3</sup> – Right (10"L x 3"W – fits ZARM10)	\$90
ZTRSL	Arm Trough <sup>4</sup> – Small Left (16" L x 4.25" W) <i>NEW</i>	\$110
ZTRSR	Arm Trough <sup>4</sup> – Small Right (16" L x 4.25" W) <i>NEW</i>	\$110
ZTRML	Arm Trough <sup>4</sup> – Medium Left (18.5" L x 5" W)	\$110
ZTRMR	Arm Trough <sup>4</sup> – Medium Right (18.5"L x 5"W)	\$110
ZTRLL	Arm Trough <sup>4</sup> – Large Left (21" L x 5.25" W) <i>NEW</i>	\$110
ZTRLR	Arm Trough <sup>4</sup> – Large Right (21" L x 5.25" W) <i>NEW</i>	\$110
ZPAML	Pivoting Arm Mount <sup>5</sup> – Left (for ZFAPLL / ZFAPLR / ZTRSL / ZTRML / ZTRLL)	\$165
ZPAMR	Pivoting Arm Mount <sup>5</sup> – Right (for ZFAPLR / ZFAPLR / ZTRSR / ZTRMR / ZTRLR)	\$165
ZASK05	Arm Support Spacer Kit <sup>6</sup> – 1/2" per side for all Adjusta-Backs (pr)	\$60
ZASK10	Arm Support Spacer Kit <sup>6</sup> – 1" per side for all Adjusta-Backs (pr)	\$115
ZPHGL	Pivoting Hand Grips (pr) for 13"L arms <sup>7</sup>	\$395
ZPHGS	Pivoting Hand Grips (pr) for 10"L arms <sup>7</sup>	\$395
ZAPSB_*	Anterior Posture Support Bar <sup>8,9</sup> * Enter frame width to complete code (14" / 16" / 18" / 20")	\$225
ZAPSB_*	Anterior Posture Support Bar <sup>8,9</sup> * Enter frame width to complete code (22" / 24")	\$245
<b>Reclining Back Arm Support Accessories</b>		
Z2PLASR	2-Point Flip-up Arm Supports (Long) (pr)	N/C
Z2PSASR	2-Point Flip-up Arm Supports (Short) (pr)	N/C
ZASKR05	Arm Support Spacers 1/2" per side (pr) for Reclining Back Support	\$60
ZASKR10	Arm Support Spacers 1" per side (pr) for Reclining Back Support	\$115
<p><sup>1</sup> Not compatible with ZPAM / ZFAPL / ZFAPS / ZTRM / ZTRL</p> <p><sup>2</sup> For chairs with same back width and frame width only</p> <p><sup>3</sup> May require ZPAM for Arm Support to flip up with Flat Arm Pad</p> <p><sup>4</sup> ZPAM is required for Arm Support to flip up with Arm Trough</p> <p><sup>5</sup> Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered.</p> <p><sup>6</sup> ZASK05 or ZASK10 is required for Standard Arm Support to flip up with lateral pad in forward position</p> <p><sup>7</sup> ZPHGL ZPHGS combined with 2 Point Arms requires a custom up charge</p> <p><sup>8</sup> Not compatible with ZFAP / ZTR</p> <p><sup>9</sup> Special length required with ZASK05 and ZASK10 (contact customer service)</p>		

**10. Foot and Leg Supports**

ZFSL	Foot Support Left – 16"-20.5"	Standard
ZFSR	Foot Support Right – 16"-20.5"	Standard
ZFET2L	Footrest Extension Tube <sup>1</sup> Left – 2" longer	\$25
ZFET2R	Footrest Extension Tube <sup>1</sup> Right – 2" longer	\$25
ZMFX8L	MFX8 Foot Support Left – 8"- 12"	\$20
ZMFX8R	MFX8 Foot Support Right – 8"- 12"	\$20
ZMFX12L	MFX12 Foot Support Left – 12"-16.5"	\$35
ZMFX12R	MFX12 Foot Supports Right – 12"-16.5"	\$35
ZADABS18L	Angle/Depth Adj ABS Footplate Left For 18"W frame	\$95
ZADABS18R	Angle/Depth Adj ABS Footplate Right For 18"W frame	\$95
ZADML	Adj Angle/Depth Molded Footplate <sup>2</sup> - Left	\$60
ZADMR	Adj Angle/Depth Molded Footplate <sup>2</sup> - Right	\$60
ZFET2AL	Footrest Extension Tube Left for ZADML – 2" longer	\$25
ZFET2AR	Footrest Extension Tube Right for ZADMR – 2" longer	\$25
ZVS_ _*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (16" / 18")	\$295
ZVS_ _*	V-Style Foot Support (1-piece) 16"-20.5" * Enter frame width to complete code (20" / 22" / 24")	\$345
ZFBAS	Flip-Back Foot Support (1-piece) 16.5"-20"	\$295
ZELRL	Elevating Leg Support <sup>3</sup> – Left (includes ZADML)	\$495
ZELRR	Elevating Leg Support <sup>3</sup> – Right (includes ZADMR)	\$495
ZRLSL	Universal Residual Limb Support – Left	\$195
ZRLSR	Universal Residual Limb Support – Right	\$195
<b>Foot Support Options for 16"W Frames</b>		
ZFS16L	Foot Support for 16"W frame Left – (16"-20.5")	N/C
ZFS16R	Foot Support for 16"W frame Right – (16"-20.5")	N/C
ZADABS16L	Angle/Depth Adj ABS Footplate Left for 16"W frame	\$98
ZADABS16R	Angle/Depth Adj ABS Footplate Right For 16"W frame	\$98
ZMFX816L	MFX8 Foot Support for 16"W Frame Left – 8"- 12"	\$20
ZMFX816R	MFX8 Foot Support for 16"W Frame Right – 8"- 12"	\$20
ZMFX1216L	MFX12 Foot Support for 16"W Frame Left – 12" - 16.5"	\$35
ZMFX1216R	MFX12 Foot Support for 16"W Frame Right – 12" - 16.5"	\$35
<sup>1</sup> Not compatible with ZMFX12/ZADM/ZVS/ZFBT		
<sup>2</sup> Not compatible with Z30016		
<sup>3</sup> ZELRL/R compatible with frame widths ≥ 18" and adjusted seat depths ≤ 18" or ≤ 21" w/ 3" longer seat frames (Z3LNGSF3)		

**11. Foot and Leg Support Accessories**

ZFSLK	Foot Support Lock (pr)	\$95
ZCS__*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (16" / 18" / 20" / 22")	\$105
ZCSIC__*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (16" / 18" / 20")	\$95
ZCSIC22	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code For 22"W frame	\$125
ZCPL	Calf Panel – Bodypoint Aeromesh <sup>1</sup> Fits frame widths (16" / 18" / 20")	\$225
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$90
ZHS2224	H-Strap – for 22" / 24"W frames – hook and loop attachment	\$110
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$60
ZNC	Neoprene Cover for Standard Footplate (pr)	\$90
ZNCA	Neoprene Cover for Angle/Depth Adjustable Footplates	\$90
ZNC16	Neoprene Covers for Narrow Molded Footplates – 16"W frame (pr)	\$90
ZNCABS16	Neoprene Cover for Narrow Molded ABS Angle/Depth Adjustable Footplates	\$90
ZNCVS	Neoprene Cover for V-Style Foot Support	\$90
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$90
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$95
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$95
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 2.5" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 ___	\$95
<sup>1</sup> Not compatible with ZMFX12 / ZFBT / ZVS		

**12. Positioning Pads** – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

<b>Head Supports</b>			
ZHSAS1	Head Support with small pad		\$495
ZHSAS2	Head Support with large pad		\$545
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment		\$75
<b>Lateral Thoracic Supports</b> – Pads on Lateral Thoracic Supports can be adjusted medially up to 3" per side. For up to 4.5" medial offset, specify ZLAEP extension plate.			
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left		\$225
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right		\$225
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left		\$225
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right		\$225
<b>Leg Adductors</b>			
ZLGAS1L	Leg Adductor w/ Small Pad – Left		\$225
ZLGAS1R	Leg Adductor w/ Small Pad – Right		\$225
ZLGAS2L	Leg Adductor w/ Large Pad – Left		\$225
ZLGAS2R	Leg Adductor w/ Large Pad – Right		\$225
<b>Hip Adductors</b>			
ZHP1L	Hip Adductor w/ Small Pad – Left		\$225
ZHP1R	Hip Adductor w/ Small Pad – Right		\$225
ZHP2L	Hip Adductor w/ Large Pad – Left		\$225
ZHP2R	Hip Adductor w/ Large Pad – Right		\$225
<b>Hip Adductors For Reclining Backs</b>			
ZHPR1L	Hip Adductor for Reclining Back w/ Small Pad – Left		\$225
ZHPR1R	Hip Adductor for Reclining Back w/ Small Pad – Right		\$225
ZHPR2L	Hip Adductor for Reclining Back w/ Large Pad – Left		\$225
ZHPR2R	Hip Adductor for Reclining Back w/ Large Pad – Right		\$225
<b>Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors</b>			
These plates allow for an additional 1.5" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors			
ZLAEPL	Lateral/Adductor Extension Plate – Left		\$28
ZLAEPR	Lateral/Adductor Extension Plate – Right		\$28
Lateral Thoracic Support	Left	Right	
Leg Adductor	Left	Right	
Hip Adductor	Left	Right	
Total Quantity ____			
<b>Medial Thigh Supports (Pommel)</b>			
ZMTSAS1	Medial Thigh Support w/ Small Pads		\$595
ZMTSAS2	Medial Thigh Support w/ Large Pads		\$595

### 13. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts				Chest Belts					
ZPBS	Pelvic Belt – Small (frame widths 14-16", 2-piece w/ side-release buckle)	\$50		ZCBS	Chest Belt <sup>2</sup> – Small (back widths 14"-16", 1-piece w/ hook and loop closure)	\$45			
ZPBM	Pelvic Belt – Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$50		ZCBM	Chest Belt <sup>2</sup> – Medium (back widths 16"-18", 1-piece w/ hook and loop closure)	\$45			
ZPBL	Pelvic Belt – Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$50		ZCBL	Chest Belt <sup>2</sup> – Large (back widths 18"-22", 1-piece w/ hook and loop closure)	\$45			
ZPBXL	Pelvic Belt – XLarge (seat frame widths 22"-24", 2-piece w/ side-release buckle)	\$65		ZCBXL	Chest Belt – XLarge <sup>2</sup> (back frame widths 22" / 24", 1-piece w/ hook and loop closure)	\$55			
ZPBICS	Pelvic Belt <sup>1</sup> – Infection Control, Small (frame widths 14"-16", 2-piece w/ side-release buckle)	\$75		ZCBICS	Chest Belt <sup>1,2</sup> – Infection Control, Small 1-piece w/ side-release buckle (back widths 14-16")	\$75			
ZPBICM	Pelvic Belt <sup>1</sup> – Infection Control, Medium (frame widths 16"-18", 2-piece w/ side-release buckle)	\$75		ZCBICM	Chest Belt <sup>1,2</sup> – Infection Control, Medium 1-piece w/ side-release buckle (back widths 16"-18")	\$75			
ZPBICL	Pelvic Belt <sup>1</sup> – Infection Control, Large (frame widths 18"-22", 2-piece w/ side-release buckle)	\$85		ZCBICL	Chest Belt <sup>1,2</sup> – Infection Control, Large 1-piece w/ side-release buckle (back widths 18"-22")	\$85			
ZPBICXL	Pelvic Belt <sup>1</sup> – Infection Control, XLarge (seat frame widths 22"-24", 2-piece w/ side-release buckle)	\$95		ZCBICXL	Chest Belt <sup>1,2</sup> – Infection Control, XLarge 1-piece w/ side-release buckle (back frame widths 22" / 24")	\$95			
ZPBBPM	Pelvic Belt – Bodypoint Aeromesh, Med	\$215		ZCBBPM	Chest Belt – Bodypoint Aeromesh, Med	\$215			
ZPBBPL	Pelvic Belt – Bodypoint Aeromesh, Large	\$215		ZCBBPL	Chest Belt – Bodypoint Aeromesh, Large	\$215			
ZPBBPXL	Pelvic Belt – Bodypoint Aeromesh, XLarge	\$175		ZCBBPXL	Chest Belt – Bodypoint Aeromesh, XLarge	\$195			
<b>Recommended sizes for Bodypoint belts</b>	Pelvic belt	ZPBBPS	N/A	ZPBBPM	16"W base frame	ZPBBPL	18"/20"W base frame	ZPBBPXL	22"/24"W base frame
	Chest belt	ZCBBPS	14" back frame	ZCBBPM	16"W back frame	ZCBBPL	18"/20"W back frame	ZCBBPXL	22"/24"W back frame
<sup>1</sup> Non-absorbent belts – recommended for infection-control environments									
<sup>2</sup> Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.									

### 14. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$67
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$125
ZSC17	Splash Curtain <sup>1</sup> for 17"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$325
ZSC20	Splash Curtain <sup>1</sup> for 20"D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1/2").	\$345
ZSCFR	Splash Curtain for 17" / 20"D (at front of aperture only, not for Front Open Seats)	\$345
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1/2"	\$345
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$345
ZBUDBH	Bard urine drainage bag holder (available for all frame widths)	\$155
<sup>1</sup> Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.		

**15. Other**

ZDCLA	Directional Caster Lock – fits on left rear caster (facilitates steering)	\$195	ZPSK1	Commode Pan Spacer Kit <sup>2</sup> – 1" lower	\$200
ATARR	Anti-Tippers <sup>1</sup> – Rear (pr)	\$245	ZPSK2	Commode Pan Spacer Kit <sup>2</sup> – 2" lower	\$215
ZCUAP4	4" Caster Upgrade (4x)	\$120	ZPSK3	Commode Pan Spacer Kit <sup>2</sup> – 3" lower	\$230
ZCUAP6	6" Caster Upgrade (4x)	\$350	ZPSK4	Commode Pan Spacer Kit <sup>2</sup> – 4" lower	\$245

<sup>1</sup> Rear Anti-Tippers are included with Reclining Back Option [Sec.7]

<sup>2</sup> Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.

**16. Credits and Omissions**

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907	Omit Foot Supports - if Foot / Leg Support is not required	-\$45
------	--	-------	------	--	-------

**Some assembly of Raz products is required.**

**17. Notes:**

--

**Specifications****Raz-AP (Attendant Propel)**

Weight Capacity	407 lbs
Overall Width	Seat width plus 4.5" (with Standard Arm Support Pads)
Overall Length (foot supports in middle position)	36.5" (with rear caster trailing), 35" (to rear canes)
Seat Depth	16" or 21" depending on seat/cushion length
Toilet Clearance Height (under seat)	
- with pan brackets	15.25"-19.25" (1" increments)
- without pan brackets	16"-20" (1" increments)
Toilet Clearance Width	16" for standard 18"W seat frame. For other widths, seat frame width minus 2"
Arm Support Height - above molded seat (1/2" increments)	
- <b>18"H Fixed Back - standard</b>	<b>8.5"</b>
- 15"H Adjusta-Back	5.5"-8.5"
- 17"H Adjusta-Back	6.5"-9.5" (1/2" increments)
- 20"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Adjusta-Back	7.5"-10.5" (1/2" increments)
- 20"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
- 24"H Reclining Adjusta-Back	7.5"-10.5" (1/2" increments)
Seat-to-Floor Height (with molded seat)	20"-24"
Product Weight (with 18"W seat)	38 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)