

## Hammock MeshPoly Sling with Head Support

### 8D2X30

- Provides support around the hips, sacral area, and under the thighs up to the shoulders with optional head support
- Divided leg supports that can be applied or removed in seated or supine positions
- Multiple loop and strap options for seated and semi-reclined positions
- Legs strap application in various positions for multi-level support
- Built-in handles for ultimate caregiver movement control




### Product Codes and Size Chart:



**Polyester Mesh (Mesh/Poly)**

Size	Part No.	A	B	C	D	W	LF3	L3	S4	HS	Max SWL	Size Trim Color
XXS	8D2730	27½" (70cm)	18" (46cm)	4½" (11cm)	4" (10cm)	25" (64cm)	19" (48cm)	23½" (60cm)	15" (38cm)	6½" (17cm)	450lbs	White
XS	8D2630	33½" (85cm)	23" (58cm)	6" (15cm)	5" (13cm)	29½" (75cm)	21½" (55cm)	23½" (60cm)	15" (38cm)	8" (20cm)	450lbs	Silver
S	8D2530	37½" (95cm)	25" (64cm)	6" (15cm)	6" (15cm)	33½" (85cm)	27" (69cm)	23½" (60cm)	15" (38cm)	10½" (27cm)	625lbs	Red
M	8D2430	42" (107cm)	27½" (70cm)	7" (18cm)	7" (18cm)	38" (97cm)	28¾" (73cm)	23½" (60cm)	15" (38cm)	13¾" (35cm)	625lbs	Yellow
L	8D2330	43" (109cm)	27½" (70cm)	8" (20cm)	8" (20cm)	41" (104cm)	32" (81cm)	23½" (60cm)	15" (38cm)	13¾" (35cm)	625lbs	Green
XL	8D2230	46½" (118cm)	29" (74cm)	10" (25cm)	8" (20cm)	46½" (118cm)	35" (89cm)	23½" (60cm)	15" (38cm)	13¾" (35cm)	1000lbs	Blue

### Specifications:



### Applications:

