

**McKESSON**

**PREVENT<sup>®</sup>G**

## Safety Winged Blood Collection Sets **GLIDE**

McKesson's Prevent G safety winged blood collection sets feature a single-use, sterile, winged needle bonded to flexible tubing with a Luer adapter. The sets are sterile and non-pyrogenic. Each blood collection set includes an easy-to-use integrated safety device that reduces accidental needle sticks. Since you don't have to change your needle insertion technique, it is easy to make the switch to Prevent G.

### FEATURES AND BENEFITS

- One-handed technique activates safety mechanism
- Compact design enhances control and flexibility
- Wing design permits improved positioning during utilization
- Flexible tubing assists in handling
- Built-in safety device
- Audible click confirms that needle is securely locked
- MFR #s 16-7BC2134 and 16-7BC2334 contain 7-inch tubing; MFR #s 16-BC2134, 16-BC2334, 16-BC2334H, 16-BC2534, and 16-BC2134H have 12-inch tubing
- MFR #s 16-BC2134H and 16-BC2334H include a pre-attached tube holder
- MFR # 16-BC2334NL comes without luer adapter
- Not made with DEHP
- Not made with natural rubber latex



### Satisfaction Guaranteed

For complete details, please visit  
[mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care<sup>™</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2024 McKesson Medical-Surgical Inc.

**MC**KESSON

**PREVENT<sup>®</sup>G**

## Safety Winged Blood Collection Sets **GLIDE**


MFR#	DESCRIPTION	TUBING	QTY
16-BC2134	Safety Winged Blood Collection Set Glide, 21G x ¾ in x 12 in	12 in Tubing	50 ea/bx, 4 bx/cs
16-BC2334	Safety Winged Blood Collection Set Glide, 23G x ¾ in x 12 in	12 in Tubing	50 ea/bx, 4 bx/cs
16-BC2534	Safety Winged Blood Collection Set Glide, 25G x ¾ in x 12 in	12 in Tubing	50 ea/bx, 4 bx/cs
16-7BC2134	Safety Winged Blood Collection Set Glide, 21G x ¾ in x 7 in	7 in Tubing	50 ea/bx, 4 bx/cs
16-7BC2334	Safety Winged Blood Collection Set Glide, 23G x ¾ in x 7 in	7 in Tubing	50 ea/bx, 4 bx/cs
16-BC2134H	Safety Winged Blood Collection Set Glide, 21G x ¾ in x 12 in	12 in Tubing	25 ea/bx, 8 bx/cs
16-BC2334H	Safety Winged Blood Collection Set Glide, 23G x ¾ in x 12 in	12 in Tubing	25 ea/bx, 8 bx/cs
16-BC2334NL	Safety Winged Blood Collection Set Glide, without Luer Adapter, 23G x ¾ in x 12 in	12 in Tubing	50 ea/bx, 4 bx/cs



**MFR# 16-BC2334**



**MFR# 16-BC2334H**



**MFR# 16-BC2334NL**

\* Need to purchase multi-sample luer adapter or syringe separately.