

## Hammock Quilted Sling with Head Support

### 5371XX

- Provides support around the hips, sacral area, and under the thighs up to the shoulders with head support
- Divided leg supports that can be applied or removed in seated or supine positions
- Multiple loop and strap options for seated and semi-reclined positions
- Legs strap application in various positions for multi-level support
- Built-in handles for ultimate caregiver movement control




### Product Codes and Size Chart:



**Quilted Mesh (Quilted)**

Size	Part No.	A	B	C	D	E	F	G	H	I	J	K	Extra Shoulder Straps	L	Max SWL	Size Trim Color
XS	537105	25" (64cm)	28" (71cm)	18" (46cm)	14" (36cm)	4½" (11cm)	4" (10cm)	5" (13cm)	19" (48cm)	10" (25cm)	19" (48cm)	16" (41cm)	537106	6" (15cm)	600lbs	Silver
S	537110	27" (69cm)	30" (76cm)	25" (64cm)	15" (38cm)	7" (18cm)	5" (13cm)	7" (18cm)	22" (56cm)	10" (25cm)	19" (48cm)	20" (51cm)	537111	6" (15cm)	600lbs	Red
M	537120	29" (74cm)	36" (91cm)	32" (81cm)	16" (41cm)	9½" (24cm)	6" (15cm)	9" (23cm)	27" (69cm)	10" (25cm)	19" (48cm)	24" (61cm)	537121	8" (20cm)	600lbs	Yellow
L	537130	31" (79cm)	42" (107cm)	39" (99cm)	17" (43cm)	12" (30cm)	7" (18cm)	11" (28cm)	32" (81cm)	10" (25cm)	19" (48cm)	28" (71cm)	537131	10" (25cm)	600lbs	Green
XL	537140	33" (84cm)	48" (122cm)	46" (117cm)	18" (46cm)	14½" (37cm)	8" (20cm)	13" (33cm)	37" (94cm)	10" (25cm)	19" (48cm)	32" (81cm)	537141	10" (25cm)	600lbs	Blue
XXL	537150	35" (89cm)	54" (137cm)	53" (135cm)	19" (48cm)	17" (43cm)	9" (23cm)	15" (38cm)	42" (107cm)	10" (25cm)	19" (48cm)	36" (91cm)	537151	12" (30cm)	600lbs	Black

### Specifications:



### Applications:

