

Features:

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- 5" front fixed wheels
- Complete with glide tips that replace rear leg rubber tips to enhance walker maneuverability on all surfaces
- Designed to meet the requirements of HCPCS code: E0143
- One (1) year limited warranty

Specifications:

	<u>Adult</u>	<u>Junior</u>
Width Inside Handgrips	18	18"
Width Inside Rear Legs	21	21"
Overall Width	22	22"
Depth Front to Rear	17	15"
Depth Folded	4	4"
Minimum Height	32	25 ½"
Maximum Height	39	32 ½"
Weight Capacity - EVENLY DISTRIBUTED	300 lb	300 lb



716270A-2	Adult	Aluminum	2/cs
716270A-4	Adult	Aluminum	4/cs
716370A-2	Junior	Aluminum	2/cs
716370A-4	Junior	Aluminum	4/cs