

PRODUCT CATALOG



AMSINO

Advancing Healthcare Worldwide®







GRAVITY ADMINISTRATION SETS

Amsino's comprehensive line of gravity administration IV sets are designed to improve safety and infection prevention for both caregivers and patients.



- Available in needle-free, pre-pierced and injection Y-sites in an array of sets to meet your clinical needs
- Available with or without pre-attached extension sets
- 10, 15 or 60 drops/mL sets available
- Includes a rotating male luer lock
- Non-pyrogenic, DEHP-free, latex-free and sterile

10 DROPS/mL




GENERAL PRIMARY SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
107201	10	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
108301	10	83 in. (211 cm)	Injection Y-site	15.6 mL	50/Case
					
108303	10	83 in. (211 cm)	Sure-Lok™ Needle-free Y-site	16 mL	50/Case
					
108304	10	89 in. (226 cm)	Non-vented, Roller Clamp (2), Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
					
109602	10	96 in. (244 cm)	(2) Injection Y-sites	17.5 mL	50/Case
					
108306	10	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	16 mL	50/Case
					
MRSA10-84	10	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	10.5 mL	50/Case
					
AC3110	10	105 in. (266 cm)	(2) Sure-Lok™ Needle-free Y-sites	20 mL	50/Case
					
MRA10E-95	10	95 in. (241 cm)	Vented/Non-vented, Pre-pierced Y-site, Slide Clamp, Pre-attached Ext Set, Slide Clamp, MicroClave® Clear Needle-free Y-site	18.5 mL	50/Case
					

10 DROPS/mL

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
108305	10	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
MRA10E-104	10	104 in. (264 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-sites, Pre-pierced Y-site, with Pre-attached Ext Set, Slide Clamp, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	13.4 mL	50/Case
					
MRSA10E-107	10	107 in. (272 cm)	Vented/Non-vented, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set	13.4 mL	50/Case
					
108308	10	83 in. (211 cm)	(2) Sure-Lok™ Needle-free Y-sites, Pre-pierced Y-site	16.6 mL	50/Case
					




FLOW CONTROLLER IV SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AFS102	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, Sure-Lok™ Needle-free Y-site	18.5 mL	50/Case
					
AFS208	10	100 in. (254 cm)	Vented/Non-vented, Flow Controller, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites	19.2 mL	50/Case
					
AFS208B	10	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, Flow Controller, (2) AMSafe® Needle-free Y-sites, Slide Clamp	19 mL	50/Case
					

10 DROPS/mL



PIGGYBACK SETS WITH BACK CHECK VALVE

- Provides you the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line








ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
1010001B	10	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, Sure-Lok™ Needle-free Y-site	19.3 mL	50/Case
					
108302	10	83 in. (211 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site, Rotating Male Luer	16.4 mL	50/Case
					
MRA102-84	10	84 in. (213 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
					
MRSA103-94	10	94 in. (239 cm)	Vented/Non-vented, Back Check Valve, Pre-pierced Y-site, (3) MicroClave® Clear Needle-free Y-sites	17.5 mL	50/Case
					
MRSA10E-110	10	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Pre-pierced Y-sites, (2) MicroClave® Clear with Pre-attached Ext Set	21.1 mL	50/Case
					

15 DROPS/mL

GENERAL PRIMARY SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AA3101	15	78 in. (198 cm)	Vented/Non-vented, Slide Clamp, Injection Y-site	15 mL	50/Case
					
158310	15	83 in. (211 cm)	Vented/Non-vented, Pre-pierced Y-site	17 mL	50/Case
					







15 DROPS/mL

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
158309	15	83 in. (211 cm)	Vented/Non-vented, Sure-Lok™ Needle-free Y-site	17 mL	50/Case
					
MRSA15-100	15	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site	12.4 mL	50/Case
					
158306	15	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	16 mL	50/Case
					
MRSA15-91	15	91 in. (231 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	11.4 mL	50/Case
					
159602	15	96 in. (244 cm)	(2) Injection Y-sites	17.7 mL	50/Case
					
MRA15-100	15	100 in. (254 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site	13 mL	50/Case
					
1510501	15	105 in. (267 cm)	Vented/Non-vented, (2) Sure-Lok™ Needle-free Y-sites	20 mL	50/Case
					
158305	15	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
					
MRSA15E-88	15	88 in. (224 cm)	Vented/Non-vented, Pinch Clamp, Pre-pierced Y-site with Pre-attached Ext Set, Pinch Clamp, (2) MicroClave® Clear Needle-free Y-sites	11.1 mL	50/Case
					
158307	15	89 in. (226 cm)	Sure-Lok™ Needle-free Y-site, Pre-pierced Y-site with Pre-attached Ext Set, Sure-Lok™ Needle-free Y-site, Slide Clamp, Pre-pierced Y-site	18.7 mL	50/Case
					
MRSA15E-120	15	120 in. (305 cm)	Vented/Non-vented, Pre-pierced Y-site, MicroClave® Clear Needle-free Y-site, Rotating Male Luer Lock, Drop-In Ext Set MicroClave® Clear Y-site	15.8 mL	50/Case
					

15 DROPS/mL

PIGGYBACK SETS WITH BACK CHECK VALVE

- Provides you the flexibility to infuse more than one bag of medication or solution
- Featuring a backcheck valve to prevent fluid exchange or backflow into the primary infusion line

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MRSA15E-96	15	96 in. (244 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Pre-attached Ext Set, Slide Clamp	12.1 mL	50/Case
					
MRA15B-105	15	105 in. (267 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites	13.23 mL	50/Case
					
AS3100	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (2) Injection Y-sites	19.7 mL	50/Case
					
AS3101	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Sure-Lok™ Needle-free Y-sites, Roller Clamp	20.8 mL	50/Case
					
AS3102	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, Slide Clamp, (2) Pre-pierced Y-sites, Roller Clamp	20.4 mL	50/Case
					
AS3110A	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Sure-Lok™ Needle-free Y-sites, Roller Clamp, Slide Clamp	22.3 mL	50/Case
					
AS3110B	15	110 in. (279 cm)	Vented/Non-vented, Back Check Valve, (3) Pre-pierced Y-sites, Slide Clamp	22 mL	50/Case
					

60 DROPS/mL

GENERAL PRIMARY SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
607201	60	72 in. (183 cm)	Injection Y-site	14 mL	50/Case
608304	60	89 in. (226 cm)	Non-vented, Roller Clamp (2), Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set with Pinch Clamp	17.5 mL	50/Case
608305	60	89 in. (226 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set, Pinch Clamp	17.5 mL	50/Case
608306	60	83 in. (211 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	16 mL	50/Case
MRSA60-84	60	84 in. (213 cm)	Vented/Non-vented, Pre-pierced Y-site, (2) MicroClave® Clear Needle-free Y-sites	10.7 mL	50/Case
MRSA60E-88	60	88 in. (224 cm)	Vented/Non-vented, Pre-pierced Y-site with Pre-attached Ext Set, (2) MicroClave® Clear Needle-free Y-sites	11.1 mL	50/Case
MRSA60-100	60	100 in. (254 cm)	Vented/Non-vented, MicroClave® Clear Needle-free Y-site, Pre-pierced Y-site	12.5 mL	50/Case

PIGGYBACK SETS WITH BACK CHECK VALVE


- Provides you the flexibility to infuse more than one bag of medication or solution
- Featuring a back check valve to prevent fluid exchange or backflow into the primary infusion line

MRA60B-92	60	92 in. (234 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites, Slide Clamp	11.4 mL	50/Case
MRSA602-100	60	100 in. (254 cm)	Vented/Non-vented, Back Check Valve, (2) MicroClave® Clear Needle-free Y-sites	12.5 mL	50/Case


60 DROPS/mL

BURETTE SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MRS60B-106	60	106 in. (269 cm)	150 mL Burette Set (w/MicroClave®, Shut-Off), (2) MicroClave® Clear Needle-free Y-sites	12.7 mL	25/Case



MRS60B-117	60	117 in. (297 cm)	24-hour 150 mL Burette Set (w/15 Micron Filter, Shut-off) with (3) MicroClave® Clear Needle-free Y-sites	13.5 mL	25/Case
-------------------	----	------------------	---	---------	---------



SELECTABLE DRIP CHAMBER AMSAFE-3® 10, 15 AND 60 DROPS/mL


AMSAfe-3® is a 3-in-1 multi-drip chamber IV set that allows providers to select a 10, 15 or 60 drops/mL setting. Designed for efficiency where inventory reduction and cost containment are essential. Now you can stock one IV set for multiple drip rates.

- Drop setting can be adjusted during medication delivery
- No need for clamps above and below chamber
- Includes male rotating luer lock
- DEHP-free, latex-free and sterile

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
D38301	10, 15 or 60	83 in. (210 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site	27.7 mL	50/Case



D38302	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok™ Needle-free Y-site	30.2 mL	50/Case
---------------	-----------------	-----------------	---	---------	---------




D38305	10, 15 or 60	91 in. (231 cm)	Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site, Back Check Valve with Pre-attached Ext Set with Pre-pierced Y-site, Slide Clamp, Sure-Lok™ Needle-free Y-site	30.5 mL	50/Case
---------------	-----------------	-----------------	---	---------	---------




SECONDARY PIGGYBACK SETS

ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AS3515	15	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	10 mL	50/Case
AS3500	20	40 in. (102 cm)	Vented/Non-vented, Roller Clamp, Hanger	9.2 mL	50/Case




Y-TYPE BLOOD SETS


ITEM CODE	DROPS/ mL	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MRSB-94	10	94 in. (239 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Slide Clamp, Roller Clamp, MicroClave® Y-site	57.2 mL	25/Case



MRSB-102	20	102 in. (259 cm)	200 Micron Filter, Roller Clamps, Hand Pump, MicroClave® Y-site	20.2 mL	25/Case
-----------------	----	------------------	---	---------	---------



MRSB-105	10	105 in. (267 cm)	170 Micron Filter, Roller Clamps, Hand Pump, Pre-pierced Y-site, Slide Clamp, Roller Clamp, MicroClave® Y-site, Pre-attached Ext Set	59.3 mL	25/Case
-----------------	----	------------------	--	---------	---------



IV EXTENSION SETS


Amsino's extension sets are available in various lengths and dimensions to accommodate clinical applications in healthcare settings. All extension sets are DEHP-free, latex-free, sterile and include a rotating male luer lock unless otherwise noted.

- Standard Extension Sets
- Microbore Extension Sets
- Smallbore Extension Sets
- Flow Controller Extension Sets
- Pressure Infusion Extension Sets


STANDARD EXTENSION SETS

Standard Extension Sets continued on page 14

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
000606	6 in. (15 cm)	Female Luer Lock, Pre-pierced Y-site, Sure-Lok™ Needle-free Y-site, Slide Clamp	2.5 mL	50/Case



STANDARD EXTENSION SETS

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
MR4T07	6 in. (15 cm)	MicroClave® Clear Needle-free Connector, T-connector	0.86 mL	50/Case
				
AE3107	7 in. (18 cm)	Removeable Sure-Lok™ Needle-free Connector, Slide Clamp	1.76 mL	50/Case
				
MR4107	7 in. (18 cm)	Removable MicroClave® Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
MR4207	7 in. (18 cm)	Bonded MicroClave® Clear Needle-free Connector, Slide Clamp	0.69 mL	50/Case
				
AE3108	8 in. (20 cm)	Removable Sure-Lok™ Needle-free Connector, Pinch Clamp	1.76 mL	100/Case
				
MR4R10	10 in. (25 cm)	Removable MicroClave® Clear Needle-free Connector, Pinch Clamp	0.98 mL	50/Case
				
003201	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Pre-pierced Y-sites	4.9 mL	50/Case
				
003202	27 in. (69 cm)	Female Luer Lock, (2) Slide Clamps, (2) Sure-Lok™ Needle-free Y-sites	5.3 mL	50/Case
				
MR4Y32	32 in. (81 cm)	(2) MicroClave® Clear Needle-free Connectors, (2) Slide Clamps	5.02 mL	50/Case
				

MICROBORE EXTENSION SETS

ITEM CODE	LENGTH	DESCRIPTION	APPROX. PRIMING VOL.	QUANTITY/ UNIT
AE0106	6 in. (15 cm)	Removeable Sure-Lok™ Needle-free Connector, Slide Clamp	0.63 mL	50/Case
				
AE1106	6 in. (15 cm)	Female Luer Lock, Slide Clamp, Fixed Male Luer Lock	0.5 mL	50/Case
		