



voyage[®]

Bon VOYAGE Adventures await

The ZIPPIE Voyage[®] combines the convenience of a quick-folding lightweight stroller frame with incredibly versatile seating and positioning, giving more options to support growing babies and their caregivers.

INTUITIVE,
EASY-TO-USE
ZIPPIE
VOYAGE
MAKES NO
COMPROMISES
WHEN IT COMES
TO WHAT'S
BEST FOR BABY



1. Canopy with Mesh Window

Keep baby shaded during the day, and open the mesh window to check in.

2. Angle Adjustable Stroller Handle

One-handed steering that adjusts to caregiver heights.

3. Quick Adjust Tilt

Effortlessly tilt forward or rearward for positioning needs.

4. Reversible Seating

Easily place the seat facing towards the caregiver, or let baby face outward to interact with their environment.

5. Transit Approved

The Voyage has been crash tested and is approved for occupied transit in a motor vehicle.

6. Carry-All Storage Basket

Store up to 15 lbs. of essentials.

7. Foot Operated Wheel Locks

Hands free! Lock wheels in place easily with your foot.

8. Airless Wheels

Flat tires are a thing of the past with airless wheels.

9. Swivel Front Wheels with Suspension and Lock-Out

Maneuver around obstacles with confidence, knowing baby is comfortable and secure. Lock-out the front wheels for smooth tracking during longer strolls.



WHEN IT COMES TO
SEATING, OUR
NUMBER ONE
PRIORITY IS
POSITIONING
YOUR CHILD FOR
HEALTHY GROWTH
AND DEVELOPMENT



Voyage Moderate Seating (VMS)

Designed for children with mild to moderate positioning needs. Offers soft, adjustable components and polka dot upholstery with colored accents.



Voyage Advanced Seating (VAS)

Designed for children with moderate to advanced positioning needs. Offers an extensive selection of seat options and adjustable components, and black stretch upholstery with colored piping.



Reversible Seat

Face the seat toward you to interact with your child, or forward so your child can explore their surroundings.

Positioning Accessories

Whitmyer® head supports, lateral trunk supports, anterior trunk supports, medial supports, and hip supports create a seating system that is entirely unique to your child.

Positioning accessories can be moved and adjusted to maintain proper support for your child as they grow.

Seating Components



T Back



Adjustable Contour Seat



Standard Planar
Seat w/hip blocks



Standard Planar Seat
w/out hip blocks

Whitmyer® Head Supports



Plush



Contour Cradle



Adjust-a-Plush

Positioning Accessories



Foot Straps



Medial Support



Butterfly Harness

Positioning Belt



Rearward (Posterior) Tilt & Recline

Improves positioning to assist with feeding, digestion, respiratory function, and visual orientation. Posterior tilt range of 0° to -30°.

Forward (Anterior) Tilt

Assists with positioning to encourage healthy development such as bringing hands to the mouth, reaching for objects, or engaging in activities. Anterior tilt range of 0° to +5°.

Height & Angle Adjustable Footrest

Grows with your child and adjusts for various needs throughout the day.

JCM Hi-Lo Base

Allows your little one to engage in activities at the table or on the floor, while being supported to meet his/her needs.



Voyage Moderate Seating Colors



Beige



Gray



Blue



Red

ACCESSORIES
TO FIT YOUR
STYLE AND
LIFESTYLE!



FUNCTIONAL ACCESSORIES

The Zippie Voyage is available with options designed to make daily tasks a breeze.



JCM Hi-Lo Base
Height Adjustable



Folding Booster Base
Height Adjustable



Upper Extremity Support Tray
Angle/depth adjustable hardware included



Grab Bar



Utility Hook



Cup Holder

Voyage Advanced Seating Colors



Purple



Pink Sorbet



Green



Navy Blue



Orange



Turquoise



Black



Burgundy



Leopard



Butterflies



Red Race Flags



Front Mount Utility Platform



Rear Accessory Platform



Battery, Vent or Utility Tray



LTV/Trilogy Ventilator Hardware
Dovetail vertical framemounted hardware



O2 Cylinder Holder



IV Support Pole



voyage®



Technical Specifications

Frame Style:

Early Intervention Device

HCPSC Code: E1232, or
E1234 (without seating)

Chair Weight: 28 lbs
(base only)

Frame Material: Aluminum

User Weight Cap.: 75 lbs

Transit Approved: Yes

Options & Accessories

Frame Type:

Tilt, Tilt & Recline

Caster Wheels: 7"

Rear Wheels:

12" Polyurethane

Moderate Seating (VMS):

One option

Advanced Seating (VAS):

Backrest Types

T-Back Cushion 1" thick

T-Back Cushion 2.5" thick

Seating Types

Standard Planar with

or without hip blocks

Contoured Wedge

Contoured Anti-Thrust

Adjustable Contour

Frame & Features

(dependent on
configuration)

Overall Width: 24"

Seat Widths:

Advanced: 6" to 13"

Moderate: 9.5" to 12"

Seat Depths:

Advanced: 5" to 13"

Moderate: 8" to 13"

Back Height:

Advanced: (Short) 7" to 13"
(Tall) 10" to 17"

Moderate: 16" to 24"

(includes head support)

Folded Height: 12"

Tilt Range: 35° (5 to -30°)