

CALCIUM ALGINATE DRESSINGS

Vitále® STRONG™ Calcium Alginate Dressings are designed to provide advanced wound care through their highly absorbent and conformable properties. Composed of strong calcium alginate fibers, these dressings are ideal for managing moderate to heavily exuding wounds. Upon contact with wound exudate or blood, the calcium alginate fibers form a hydrophilic gel that maintains a moist wound environment, which is crucial for promoting faster healing and reducing the risk of maceration to the surrounding skin. The dressings are highly absorbent, extending wear time, and they can be removed in one piece during dressing changes without disturbing granulating tissue. This makes Vitále® STRONG™ Calcium Alginate Dressings an excellent choice for both superficial and deep wounds, including cavity or tunneled wounds.



USAGE

Type: Primary

Stages: II-IV

Drainage/Exudate: Moderate to Heavy

Wound Types:

Pressure Sores

Venous Ulcers

Arterial Ulcers

Diabetic Ulcers

Grafts & Skin Donor Sites

Surgical/Post-Op

Traumatic Wounds

FEATURES & BENEFITS

- **Highly Absorbent:** The dressings efficiently absorb exudate, forming a gel that maintains an optimal moist wound environment for healing.
- **Non-Adherent:** The gel-like covering prevents the dressing from sticking to the wound bed, ensuring painless removal and minimizing trauma.
- **Reduces Risk of Maceration:** The highly absorbent nature of the dressing helps protect the surrounding skin by reducing the risk of maceration.
- **Conformable and Flexible:** The soft and flexible design allows the dressing to conform to the wound's shape, making it ideal for use in cavity or tunneled wounds.
- **Versatile Application:** Suitable for a wide range of wounds, including pressure sores, arterial ulcers, venous ulcers, diabetic ulcers, and post-operative wounds.

Vitále® STRONG™ Calcium Alginate Dressings are trusted by healthcare professionals for their advanced wound management capabilities. The combination of high absorbency, non-adherent properties, and flexibility ensures that wounds are kept in an optimal moist environment, promoting faster healing and reducing the risk of maceration. These dressings are versatile, easy to apply, and suitable for a wide range of wound types, making them a critical tool in both clinical and home care settings. By choosing Vitále® STRONG™ Calcium Alginate Dressings, healthcare providers can deliver superior wound care that enhances patient outcomes and satisfaction.

ORDERING INFORMATION

CATALOG NUMBER	DESCRIPTION	SIZE	UNITS/BOX	HCPC CODE
20650	Vitále® STRONG™ Calcium Alginate Rope	3/4" x 12" (2 cm x 30 cm)	5	A6199
20652	Vitále® STRONG™ Calcium Alginate Dressing	2" x 2" (5 cm x 5 cm)	10	A6196
20654	Vitále® STRONG™ Calcium Alginate Dressing	4" x 4" (10 cm x 10 cm)	10	A6196
20655	Vitále® STRONG™ Calcium Alginate Dressing	4" x 5" (10 cm x 12.5 cm)	10	A6197
20656	Vitále® STRONG™ Calcium Alginate Dressing	6" x 6" (15 cm x 15 cm)	5	A6197
20658	Vitále® STRONG™ Calcium Alginate Dressing	4" x 8" (10 cm x 20 cm)	5	A6197