



**A. Model 7011**  
Shown with option:  
32 – Pillow



**B. Model 7031**  
Shown with options:  
05 – Headrest  
27 – Paper Dispenser and Cutter

## Common Features

RECOVERY COUCHES —  
MODELS 7011, 7031, 7155, 7905:



- Wood has no added urea-formaldehyde for safer indoor air quality.
- PVC-Free construction.
- 2" High-Density urethane foam top.

**COLORS:**

Upholstery – 5 PVC-Free Urethane colors.  
(FOLD OUT FRONT COVER)

## OPTIONS

FOR MODELS 7011, 7031, 7155, 7905

- \*04 – Table Safety Strap
- \*05 – Headrest, adj. (7011, 7031 only)
- \*07 – Nose Cutout
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small

\*FACTORY INSTALLED

SEE PAGE 4

### A. Model 7011 Econo-Line™ Recovery Couch



#### SPECIAL FEATURES:

- Sturdy end panel legs and apron.
- This stylish and functional couch is ideal for any school or health facility.
- Folkstone Gray laminate finish.
- Shipped ready-to-assemble.

MODEL	LENGTH	WIDTH	HEIGHT
7011	72"	24"	22"

### B. Model 7031 Econo-Line™ Recovery Storage Couch



- Cabinet with (2) sliding laminate doors on front and back side. Provides convenient storage space for all blankets, gowns, and supplies.
- Folkstone Gray laminate finish.
- Shipped ready-to-assemble.

MODEL	LENGTH	WIDTH	HEIGHT
7031	72"	24"	22"

### C. Model 7155 Fleetwood Couch



#### FEATURES:

- (10) position, 0-45° adjustable Headrest.
- Fusion Maple laminate legs.

MODEL	LENGTH	WIDTH	HEIGHT
7155	72"	27"	18"

### D. Model 7905 Inwood Couch



#### FEATURES:

- 3" thick urethane foam vinyl top.
- (10) position, 0-45° adjustable Headrest.
- Fusion Maple laminate legs.

MODEL	LENGTH	WIDTH	HEIGHT
7905	72"	27"	18"



**C. Model 7155**



**D. Model 7905**