

# Protekt® BP

## Upper Arm Blood Pressure Monitor



HCPCS  
A4670



### FEATURES AND BENEFITS:

- Accurate result detected by oscillograph.
- Easy one-touch operation.
- Large, easy to read white backlit LCD screen.
- Displays pulse, systolic and diastolic pressure.
- Irregular heartbeat detection (IHD).
- WHO classification indicator.
- Dual language broadcast, English & Spanish voice feature.
- 80x3 users - memory storage (240 total memory recall).
- Average of the 3 most recent measurements.
- Date & time display.
- Digital error messages.
- Movement error detection.
- Cuff wearing detection.
- Automatically powers off to help extend battery life.
- USB type-C port for power charge.
- Unit dimensions: 4.65"(L) x 3.74"(W) x 2.6"(H)
- LCD screen size: 2.83"(L) x 2.2"(W)
- Arm cuff circumference: Fits 8.7" to 16.5"
- Storage bag included.
- Works with 4 "AA" batteries (included).
- Warranty: 2 Year limited.

PMDBPA

Protekt® BP Blood Pressure Monitor

