

# New generation of Molift mobile lifts


Caring for others requires equipment you can depend on, and for decades, Molift mobile lifts have been recognized globally for their unwavering performance. To ensure we continue to meet your evolving needs, we've listened, innovated, and are proud to introduce our upgraded Molift floor lift series. Experience a new level of functionality and comfort, as these lifts are now equipped with advanced features designed to optimize your experience and provide lasting value.



**New battery with increased capacity**



**Significant increase of product lifetime**



**Bluetooth connection to Molift Service Tool**



**Turbo boost when unloaded**

## Upgraded lithium-ion battery technology

With modern lithium-ion battery technology the new generation of Molift lifts offers more than double the capacity and over three times the life expectancy compared to previous models.

Additionally, they now come with a five-year warranty, an improvement from the current one-year warranty. With a life expectancy of 1,700 charging cycles, compared to the old version's 500 cycles, the battery is built to last. The new battery is also lighter in weight (2 lbs.) and the same model is used for all new Molift lifts.

## Increase of product lifetime

The upgraded lifts are even more durable and sustainable than the previous versions, with the expected product lifetime increased from 10 to 15 years.

## Bluetooth connection to Molift Service Tool

The new lifts enable wireless Bluetooth access to Molift Service Tool to conduct service and access data on hoist usage, lifting patterns etc.

## Turbo boost function

When used without load, the lift goes **30% faster**.

## New hand control

More robust and washable (IPX6) hand control with LED indicator for battery level, overload and service.



**molift®**

## New article numbers

Below, you'll find a comprehensive list of all new part numbers, along with their corresponding legacy products.

Product Name	Description	New Part Number	Legacy Part Number
Molift Smart 150	Standard Molift Smart 150	N1020	N10200
Molift Mover 180	Standard model, 3-inch (75mm) wheels	N1801	N18000
	Standard model, 4-inch (100mm) wheels	N1802	N18003
Molift Mover 205	Standard model	N1401	N14000
	Low Base	N1402	N14050
	Standard model, includes support arms	N1403	N14150
Molift mover 300	Standard model	N1501	N15050
	Low Base	N1502	N15051
	Standard model, includes support arms	N1503	N15150
	Low Base, includes support arms	N1504	N15151