



# Manatee

Company:	<input type="text"/>	Ship To:	<input type="text"/>
Account Number:	<input type="text"/>	Date:	<input type="text"/>
<input type="checkbox"/> Quote: <input type="checkbox"/> PO Number:	<input type="text"/>		
ATP:	<input type="text"/>	City:	<input type="text"/>
		State:	<input type="text"/>
Marked For:	<input type="text"/>	Zip:	<input type="text"/>
		Phone:	<input type="text"/>
Associated Quote/Order#:	<input type="text"/>	Email/Fax:	<input type="text"/>

Effective: 1/5/2026

## Manatee

	mm (inch)	kg (lbs)		
		0		1
		2		3
Effective seat width	330 mm (13")	330 mm (13")	330 mm (13")	400 mm (15¾")
Effective seat depth	220-250 mm (8¾-9¾")	300-370 mm (11¾-14½")	390-460 mm (15¼-18")	400-470 mm (15¾-18½")
Backrest height	430 mm (17")	430 mm (17")	580 mm (22¾")	640 mm (25¼")
Lower leg length	280 mm (11")	280 mm (11")	350 mm (13¾")	350 mm (13¾")
Seat surface height at front edge	140-430 mm (5½-17")	140-430 mm (5½-17")	140-430 mm (5½-17")	140-430 mm (5½-17")
Overall width	420 mm (16½")	420 mm (16½")	420 mm (16½")	490 mm (19¼")
Folded length	1030-1100 mm (40½-43¾")	1100-1170 mm (43-46")	1360-1430 mm (53½-56¾")	1440-1510 mm (56¾-59½")
Folded height	130 mm (5")	130 mm (5")	130 mm (5")	130 mm (5")
Mass device	5.5 kg (12.1 lb)	5.5 kg (12.1 lb)	6 kg (13.2 lb)	7 kg (15.4 lb)
Maximum user mass	50 kg (110.2 lb)	50 kg (110.2 lb)	75 kg (165.3 lb)	75 kg (165.3 lb)
Maximum rated load	57.5 kg (126.7 lb)	57.5 kg (126.7 lb)	83 kg (182.9 lb)	84 kg (185.1 lb)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1259.2 Pricing and specifications subject to change without notice.



## Manatee Seat and Frame Options

### Size and Color

*Shower seat only*



#### Blue

- ☐ Size 0, Blue (882005) \$780
- ☐ Size 1, Blue (882015) \$780
- ☐ Size 2, Blue (882025) \$780
- ☐ Size 3, Blue (882035) \$852

#### Red

- ☐ Size 1, Red (882014) \$780
- ☐ Size 2, Red (882024) \$780
- ☐ Size 3, Red (882034) \$852

### Manatee Frame Options



#### Rolling Shower Frame

- ☐ Rolling Shower Frame (882100) \$2499



#### Tub Slider Frame

- ☐ Tub Slider without Manatee W/ Casters, Size 1 and 2 (MC6000PED) \$4100
- ☐ Tub Slider without Manatee W/ Casters, Size 3 (MC6200PED) \$5100



#### Basic Manatee Frame

- ☐ Basic Manatee Frame (882101) \$390

## Accessories and Straps

### Accessories



#### Leg Straps

- ☐ Leg Straps, Size 0,1 and 2, Blue (883055) \$86
- ☐ Leg Straps, Size 0,1 and 2, Red (883054) \$86
- ☐ Leg Straps, Size 3, Blue (883105) \$109
- ☐ Leg Straps, Size 3, Red (883104) \$109



#### Head/ Lateral Support

- ☐ Head Lateral Support, Size 0,1 and 2, Blue (883015) \$123
- ☐ Head Lateral Support, Size 0,1 and 2, Red (883014) \$123
- ☐ Head Lateral Support, Size 3, Blue (883065) \$123
- ☐ Head Lateral Support, Size 3, Red (883064) \$123



#### Pommel

- ☐ Pommel, Size 0,1,2, Blue (883045) \$109
- ☐ Pommel, Size 0,1,2, Red (883044) \$109
- ☐ Pommel, Size 3, Blue (883095) \$109
- ☐ Pommel, Size 3, Red (883094) \$109



#### Torso Belt

- ☐ Torso Belt, Size 0,1 and 2, Blue (883025) \$86
- ☐ Torso Belt, Size 0,1 and 2, Red (883024) \$86
- ☐ Torso Belt, Size 3, Blue (883075) \$86
- ☐ Torso Belt, Size 3, Red (883074) \$86