

# Delta Terry-Net™

LATEX  
FREE  
FORMULA

## Adhesive Terry Cloth/Foam Padding

Delta Terry-Net™ Adhesive Cloth/Foam Padding provides added protection and cushioning by combining a layer of soft terry cloth bonded to 1/8" foam. The adhesive backing makes attaching the padding simple and easy.

### Features & Benefits

- Open cell foam
- Adhesive back
- Comfortable terry cloth cushion

#### Delta Terry-Net™ Thumb Spica Adhesive Terry Cloth/Foam Padding

| Ref. No. | Size          | Unit of Measure |
|----------|---------------|-----------------|
| 55012    | 23 in x 39 in | 1 roll / case   |



# Delta Terry-Net™

LATEX  
FREE  
FORMULA

## Adhesive Felt Padding

Delta Terry-Net™ Adhesive Felt Padding is an adhesive orthopedic felt, used for added patient comfort. Its adhesive backing makes attaching the padding simple and easy. It may also be used as a cast edge padding or to provide protection for bony prominences.

### Features & Benefits

- Soft orthopedic felt
- Adhesive back

#### Delta Terry-Net™ Adhesive Felt Padding

| Ref. No. | Size          | Unit of Measure |
|----------|---------------|-----------------|
| 54012    | 23 in x 39 in | 1 roll / case   |

