

Rose Comfort Max tilt and recline chair with casters

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- The DRIVE Chair with Caster Full Tilt and Recliner offers similar features but at a competitive, entry level price point. Perfect for the LTC facilities and rental companies looking for value product. When stacked against our competition, it is the only value priced product on the market with a 400lb weight capacity. The DRIVE Chair with Caster Full Tilt and Recliner will offer exceptional positioning and comfort over traditional wheel and geriatric chairs.
- 38 degrees of tilt- lever on the back left positioned easily for the caregiver to adjust
- Recline up to 90 degrees (trendelenburg). Recline lever is Red and located on the back Right
- Seat Width 20" and seat depth 20"
- Seat height adjusts from 17"-22" in 1" increments
- Leg rest: Elevating, length adjustable, lever in the back or side positioned easily for caregiver to adjust
- Footrest: Removable single flip up footrest with angle adjustable foot plate, need also secure so it does not get lost
- Wings: Adjustable Removable for Upper Lateral Support
- Casters: (2) Front 5 total lock (12.5 cm) Rear (1) 5 Directional locking & (1) 5 total lock caster (12.5 cm)
- Standard Padding: Headrest cushion, shoulder bolsters, lower side pads, upper side pads, seat and back pad with cover, calf pad, sole pad

PRODUCTS

Description

CLN400-20	Chair with Caster Full Tilt and Recliner
CLN400-HPP	Huntington's Padding Package CLN400-20

SPECIFICATIONS

Description

Warranty	Lifetime (10 years on frame)
Weight Capacity	400 lbs
Color	Blue
Dimensions	Arm Pad: 14" (L) x 3" (W) Leg Pad: 19" (L) x 12.5" (H) Foot Pad: 12" (W) x 22" (H)
Overall Width	25.8"
Recline	Up to 90 degrees (trendelenburg)
Tilt Angle Range	Up to 38 degrees
Back Height	33"

ACCESSORIES

Description

CLN400-HPP	Huntington's Padding Package CLN400-20, EA
------------	--