

## In-Situ High Easy Superfine Plus Clips

### M921530

- Designed to provide full body and head support on transfer
- Adjustable Integral head support
- Easy to apply with low-friction fabric to leg areas for ease of application
- Soft binding as standard
- Universal clip fixing
- Once the resident has been transferred, this sling can be left in-situ, subject to an appropriate risk assessment
- Wash temp 95°C / 203°F
- Tumble dry low heat (max 50°C)



### Product Codes and Specifications:



Silva-Superfine PLUS®		A	B	C	D
XS	M921530-XS	30 - 34" (76 - 86cm)	23¾ - 26" (60 - 66cm)	26¾" (68cm)	27¼" (69cm)
S	M921530-S	34 - 38" (86 - 96cm)	26 - 29¼" (66 - 74cm)	30" (76cm)	30¼" (77cm)
M	M921530-M	38 - 42½" (96 - 108cm)	29¼ - 31½" (74 - 80cm)	34" (86cm)	33" (84cm)
L	M921530-L	42½ - 51¼" (108 - 130cm)	31½ - 34¾" (80 - 88cm)	38½" (98cm)	36¾" (93cm)
XL	M921530-XL	51¼ - 56¾" (130 - 144cm)	34¾ - 39½" (88 - 100cm)	40¼" (102cm)	41" (104cm)
XXL	M921530-XXL	56¾ - 63" (144 - 160cm)	37¾ - 41" (96 - 104cm)	47¼" (120cm)	42½" (108cm)

